

Magnetic recumbent exercise bike for home fitness. This exercise bike is ideal for lower limb muscular recovery as the tension is manually adjustable. The most important details of the workout are constantly shown on the display. The user seats in a more comfortable position compared to usual exercise bikes.



Handgrip

TECHNICAL SPECIFICATIONS

Product code	1811030
Resistance System	Magnetic
Tension Levels	8 levels manual tension
Flywheel	9 kg
Display	1 window lcd
Functions	Time / speed / distance / calories / pulse
Programs	No
Heart Rate Monitoring	Hand pulse
Cardio belt	No
Handlebar	Ergonomic, antislip
Seat	High absorption, horizontal adjustable, breathable backrest
Crank arms	3 pcs
Maximum User Weight	135 kg
Net Weight	43 kg
Gross Weight	49 kg
Set-Up Sizes	152 x 60,5 x 125,5 cm (LxWxH)
Package Dimensions	150 x 28 x 69 cm (LxWxH)
Transportation Wheels	Yes
Certifications	CE, ROHS, EN957

FISIOTECH SRL

Sede Operativa e Amministrativa: Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

Sede legale : Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail: info@fisiotech.com - pec: fisiotech@emailsec.it

P.I. e C.F. 03765100262 - Rea TV 296956