


Cut-out wedge cushion.  
Polyurethane foam padding density: 30 kg/m<sup>3</sup>, Italian fire  
resistance rating: 1 IM (the highest). Washable leatherette  
upholstery, Italian fire resistance rating: 1 IM.  
Colour blue.



## MODELS – DIMENSIONS TECHNICAL SPECIFICATIONS

**L**  
40 cm

**P**  
34 cm

**h**  
9 cm

**Code**  
1415407

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956

## TECHNICAL SPECIFICATIONS – PADDING

|                        |                      |
|------------------------|----------------------|
| Name of material       | Polyurethane foam    |
| Reaction to fire class | 1 IM (highest)       |
| Density                | 30 kg/m <sup>3</sup> |
| Height                 | 40 mm                |



## TECHNICAL SPECIFICATIONS – UPHOLSTERY

| FEATURE                                                  | TEST METHOD                                                                                                                                  | UNIT          | LIMIT VALUE                                            |
|----------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------------------------------------------------------|
| Tensile strength                                         | UNI EN ISO 1421<br>Method 1                                                                                                                  | N             | Lengthwise ≥ 210<br>Crosswise ≥ 90                     |
| Elongation strength                                      | UNI EN ISO 1421<br>Method 1                                                                                                                  | %             | Lengthwise ≥ 85<br>Crosswise ≥ 140                     |
| Determination of tear force                              | EN ISO 13937-2: 2000                                                                                                                         | daN           | Lengthwise ≥ 2,5<br>Crosswise ≥ 2,2                    |
| Colour fastness to crocking                              | UNI EN ISO 105 X12: 2003<br>UNI EN ISO 105 D02: 1998                                                                                         | grey<br>scale | Dry 4<br>Bleach 4/5<br>Ethyl alcohol 3/4               |
| Determination of abrasion resistance: Martindale abrader | UNI EN ISO 5470-2: 2004                                                                                                                      | cycles        | Method 1 (direct) 65.000<br>Method 2 (reverse) 210.000 |
| Reaction ignition<br>Fireproof class                     | UNI 9175: 1987<br>UNI 9175/FA1: 1994                                                                                                         |               | CLASSE 1.IM                                            |
| Free and hydrolyzed formaldehyde                         | UNI EN ISO 14184-1: 2011                                                                                                                     | mg/Kg         | < 3                                                    |
| Phtalates                                                | UNI EN ISO 14389: 2014                                                                                                                       | %             | < 0,06                                                 |
| Arylamines                                               | UNI EN ISO 14362-1: 2012                                                                                                                     | mg/Kg         | Not detectable                                         |
| Cleaning                                                 | Clean with damp soapy cloth and rinse well with clean water. Do not use solvents, bleaches, synthetic detergents, polishes or aerosol spray. |               |                                                        |
| Specific features                                        | Antibacterial Fingicide                                                                                                                      |               |                                                        |
| Composition                                              | Plasticised PVC 88% - Polyester 12%                                                                                                          |               |                                                        |

**FIRE RESISTANCE ACCORDING TO:**  
IMO - 2010 FTP CODE Part 8  
"Test for upholstered furniture"  
par 7.2 SMOULDERING CIGARETTE TEST  
par 7.3 PROPANE FLAME TEST  
EN 1021-1: 2006 smouldering cigarette test  
EN 1021-2: 2006 match flame equivalent test

### FISIOTECH SRL

**Operating headquarters:** Via Toniolo, 44/E-F - 31028 Vazzola (TV) - ITALY

**Registered office:** Via Valdoni, 19 - 31013 Cimavilla di Codognè (TV) - ITALY

Tel. +39 0438470342 - Fax +39 0438478762 - mail : [info@fisiotech.com](mailto:info@fisiotech.com) - pec : [fisiotech@emailsec.it](mailto:fisiotech@emailsec.it)

VAT Number IT 03765100262 - Rea TV 296956