FISIOTECH

TECHNOLOGIES AT THE SERVICE OF HEALTH

since
2001

SINCE 2001 ... together for more than twenty years

Fisiotech was established in 2001 with the idea of manufacturing quality rehabilitation equipment. Helping, through new ideas, professionals in performing their precious work of prevention and care of the person: this has always been our objective.

This simple, fundamental idea is the basis of our payoff:

TECHNOLOGIES AT THE SERVICE OF HEALTH

Our operating headquarters is also the place where our ideas are born, where we create the projects, develop the prototypes, and lastly build our products.

The sales, administration and technical departments are available to develop and fabricate custom products based on the customer's specifications.

The internal production departments, from carpentry to upholstery, from assembly to warehouse allow us to respond promptly to the main requests of our customers.

We entrust all design, production, and assembly steps to our qualified staff, who carry out rigorous checks in the different manufacturing steps.

The application of a certified quality system guarantees the quality of the products manufactured by FISIOTECH. Specifically, our company holds the certifications according to the International standards UNI EN ISO 9001 and UNI EN ISO 13485 and complies with the European Regulation EU 745/2017 (MDR) on the safety of medical devices.


Fisiotech would like to thank all those who contributed to the realization of this catalogue, for their work, their advice, and their support.



THE NEW PARAMETRIC CODE

HOW TO CREATE IT ACCURATELY TO ORDER COUCHES

The will to offer more performances and customization has led us to add new options to our couches.
This sometimes generates difficulties when configuring the product as it was ordered by the customer.
The parametric code helps the customer to compose in an easy and univocal way a code for the chosen couch containing all the desired options.



Sometimes this section, too, could contain preset values.
This means that the configuration cannot be modified.

EXAMPLE:
How to order a Giove 5 RAL 7035 LARGE, with all-around foot bar switch and 2 actuators, comfort padding, Burgundy colour.

PARAMETRIC CODE

1	1	1	8	E	3	1	C2	C	B2
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PEDAL 2 ACTUATORS	ALL AROUND 2 ACTUATORS	ALUM. ALL AROUND 2 ACTUATORS
P2	C2	A2
8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

YOUR HYPOTHERTICAL CHOICE

BIT 6 RAL 7035 painting = 3

BIT 7 couch width 80 cm = 1

BIT 8 all-around foot bar switch and 2 actuators movement = C2

BIT 9 comfort padding = C

BIT 10 upholstery colour Burgundy = B2 as shown in chart.

INDEX

COUCHES

Top Series....5
Professional Series....21
Simple Series....33
Bobath Series....41
Traction Series....53
Tilt Series....61
Postural Series....71
Set-Close Series....79
Gynecology Series....89
Easy Series....95
Upholstery....103
Top Series Variations and Accessories....104
Professional Series Variations and Accessories....105
Simple Series Variations and Accessories....106
Bobath Series Variations and Accessories....107
Traction Series Accessories....108
Tilt Series Variations and Accessories....109
Postural Series Variations and Accessories....110
Set-Close Series Accessories....112
Gynecology Series Variations and Accessories....113
Easy Series Variations and Accessories....114
General Accessories....115

REHABILITATION

Parallel Bars and Accessories....119
Staircases and Accessories....124
Rocker Boards....126
Movin'step – Stability Wheel – Core Balance....127
Disco Sport – Disco Dome – Accessories....128
Artemide Limb Rehabilitation Armchair – Pedal Exerciser....129
Shoulder Wheels....130
Hand Ladder - Spirals – Worm Screw – Ergotherapy Board....131
Hand Grip – Foam Balls –
Wooden Spheres and Disks – Gas Pedal Footboard....132
Fisiotek and Accessories....133
Pulley Therapy Cage and Accessories....136
Pulleys – Pulley with Wall Bar and Accessories....139

GYMNASTICS

Gymnic Classic and Plus....143
Therasensory - Physio Roll....144
Foam Balls – Softgym – Medicine Balls....145
Rhythmic Balls – Sensy-Balls....146
Reflex Massage Balls – Thera Freeballs Hand – Activ Roll....147
Accessories for Balls....148
Wall Bars and Accessories....149
Mirrors and Accessories....152
Mattresses and Accessories....153
Mats....154
Airex Mats and Accessories....156
Balance Pad – Balance Beam Airex....162
Roll and Half-roll Cushions....163
Cube and Parallelepiped Cushions....164
Wedge Cushions....165
Trapezoid Cushions - Movin'sit....166
Disco'sit and Accessories....167
Hoops – Sticks – Cones with holes....168
Clubs – Wooden Clubs and Accessories....169
Elastics and Accessories....170
Wrist and Ankle Weights....173
Dumbbells and Accessories....174
Fit Ball Roller – Foam Roll and Accessories....175
Skipping Ropes – Mini Trampoline –
Baumann Blocks and Supports – Magic Circle....176
Wooden Trolley and Psychomotor Kit....177

DIAGNOSTICS

Podoscope – Steps for stress test –
Plumb Bob – X-Ray Film Viewers....181
Eye Charts – Scales – Height Rod....182

FURNITURE

Glass Door Cabinets – Desk –
Revolving and Upholstered Chair....185
Coat Rack – Screen....186
Revolving Stools – Upholstered Benches....187
Wooden Stools – Wooden Chairs – Medical Trolley....188
Stainless Steel Trolley – Couch Stairs....189
Changing Room Benches....190
Wall-mounted Clothes Racks – Footrests....191

ORTHOPAEDIC AIDS

Crutches – Walkers....195
Walkers – Rollator....196

HOME FITNESS


Easy Access Exercise Bike....199
Electromagnetic Exercise Bike....200
Recumbent Exercise Bike....201
Treadmill....202
Stepper....203
Rowing Machine....204

Terms of Sale....206

BREATHING HOLE

STANDARD


S



STANDARD

COMFORT

C



COMFORT

UPHOLSTERY COLOURS

Beige P1



Light grey G2



Taupe P2



Dark grey G3



Yellow G1



Mustard G4



Orange A1



Red R1



Burgundy B2



Violet V4



Dark blue B3



Blue B1



Dark brown M1



Earth M2



Green V1



Forest green V2



Leaf green V3




Black N1

FISIOTECH


TOP

FISIOTECH






TOP




PROFESSIONAL EXAMINATION
AND TREATMENT COUCHES



TOP SERIES

Couches with a unique frame that allows height adjustment without lateral offset.





GIOVE 1 | TOP

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1118A30C1SA1 + Accessories Code 1218020







Standard padding



Comfort padding

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	190 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	A					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY


10 — See colours and materials in the chart on the back cover flap.


GIOVE 2 | TOP

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1118B20C1SB2 + Accessories Code 1218020







Aluminium all around foot bar switch



Comfort padding

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	B					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

GIOVE 3 | TOP

3-SECTION ELECTRICAL OR HYDRAULIC COUCH




Model shown Code 1118C20C1SV1 + Accessories Code 1218020





Hydraulic control



TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Middle section adjustment	with gas spring or electrically +20°
Footrest adjustment	with gas spring +85°
Trendelenburg position	+20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	C					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PEDAL 2 ACTUATORS	ALL AROUND 2 ACTUATORS	ALUM. ALL AROUND 2 ACTUATORS
P2	C2	A2
8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

GIOVE 4 | TOP

4-SECTION ELECTRICAL OR HYDRAULIC COUCH




Model shown Code 1118D30C1SG2 + Accessories Code 1218020




Manual armrest adjustment



All around foot bar switch





TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Armrests on headrest adjustment	manual
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	D					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10 — See colours and materials in the chart on the back cover flap.

GIOVE 4 PLUS | TOP

4-SECTION ELECTRICAL OR HYDRAULIC COUCH




Model shown Code 1118H30C1SG2 + Accessories Code 1218020





Manual armrest adjustment



TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +80° -20°
Armrests on headrest adjustment	manual
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	H		0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016 — 2 — 6

RAL 7035 — 3 — 6

WIDTH

65 cm — 0 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

ALUM. ALL AROUND 1 ACTUATOR — A1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PADDING

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

10 — See colours and materials in the chart on the back cover flap.

GIOVE 5 | TOP

5-SECTION ELECTRICAL OR HYDRAULIC COUCH




Model shown Code 1118E20C2SB2 + Accessories Code 1218020



Trendelenburg position



Comfort padding



TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Armrests on headrest adjustment	manual
Middle section adjustment	with gas spring or electrically +20°
Footrest adjustment	with gas spring +85°
Trendelenburg position	+20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	E					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016 — 2 — 6

RAL 7035 — 3 — 6

WIDTH

65 cm — 0 — 7

80 cm — 1 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

ALUM. ALL AROUND 1 ACTUATOR — A1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PEDAL 2 ACTUATORS — P2 — 8

ALL AROUND 2 ACTUATORS — C2 — 8

ALUM. ALL AROUND 2 ACTUATORS — A2 — 8

PADDING

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY


10


See colours and materials in the chart on the back cover flap.

GIOVE 6 | TOP

6-SECTION ELECTRICAL OR HYDRAULIC COUCH







Model shown Code 1118F20A1SA1 + Accessories Code 1218020



Manual armrest adjustment



Aluminium all around foot bar switch

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Armrests on headrest adjustment	manual
Removable armrests rotating by	180°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Double-ply paper roll	1210060
Safety magnet	1218030

ACCESSORIES	CODE
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	F		0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm
0
7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

GIOVE 7 | TOP

7-SECTION ELECTRICAL COUCH




Model shown Code 1118G30C2SN1 + Accessories Code 1218020





Sitting position



Trendelenburg position







TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 49-95 cm
Length	200 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Armrests on headrest adjustment	manual
Middle section adjustment	electrically +20°
Middle section armrests	foldable
Footrest adjustment	with gas spring +85°
Trendelenburg position	+20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE	ACCESSORIES	CODE
Paper roll holder	1218010	Safety magnet	1218030
Four castors	1218020	Safety double click	1218040
Terry towel couch cover 65 cm	1210040	Swiss plug	1218050
Double-ply paper roll	1210060	UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	G		0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016 — 2 — 6

RAL 7035 — 3 — 6

WIDTH

65 cm — 0 — 7

MOVEMENT

PEDAL 2 ACTUATORS — P2 — 8

ALL AROUND 2 ACTUATORS — C2 — 8

ALUM. ALL AROUND 2 ACTUATORS — A2 — 8

PADDING

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

10


See colours and materials in the chart on the back cover flap.


ZEFIRO 2 | TOP

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1118M30C1SR1 + Accessories Code 1218020







Trendelenburg position



All around foot bar switch

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Trendelenburg position	with gas spring or electrically +20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

PARAMETRIC CODE

1	1	1	8	M					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PEDAL 2 ACTUATORS	ALL AROUND 2 ACTUATORS	ALUM. ALL AROUND 2 ACTUATORS
P2	C2	A2
8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10


See colours and materials in the chart on the back cover flap.


ZEFIRO 4 | TOP

4-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1118N20C1SA1 + Accessories Code 1218020







Trendelenburg position



Hydraulic control

PARAMETRIC CODE

1	1	1	8	N					
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016	RAL 7035
2	3
6	6

WIDTH

65 cm	80 cm
0	1
7	7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PEDAL 2 ACTUATORS	ALL AROUND 2 ACTUATORS	ALUM. ALL AROUND 2 ACTUATORS
P2	C2	A2
8	8	8

PADDING

STANDARD	COMFORT
S	C
9	9

UPHOLSTERY

10 — See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Armrests on headrest adjustment	manual
Trendelenburg position	with gas spring or electrically +20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.


ZEFIRO 6 | TOP

6-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1118P30C1SB1 + Accessories Code 1218020





PARAMETRIC CODE

1	1	1	8	P		0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7016 — 2 — 6

RAL 7035 — 3 — 6

WIDTH

65 cm — 0 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

ALUM. ALL AROUND 1 ACTUATOR — A1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PADDING

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

10 — See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 49-95 cm
Length	200 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +70° -40°
Armrests on headrest adjustment	manual
Middle section	fixed
Upper footrest adjustment	with gas spring +50°
Lower footrest adjustment	with gas spring +15° -52°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	90 kg

ACCESSORIES	CODE
Paper roll holder	1218010
Four castors	1218020
Terry towel couch cover 65 cm	1210040
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1218030
Safety double click	1218040
Swiss plug	1218050
UK 3-pin plug	1218051

ACCESSORIES AND VARIANTS see descriptions and pictures on page 104.

FISIOTECH

TECHNOLOGIES AT THE SERVICE OF HEALTH



NEW PROFESSIONAL

BREATHING HOLE

STANDARD

S





STANDARD

COMFORT

C





COMFORT

UPHOLSTERY COLOURS

Beige P1

Light grey G2

Taupe P2

Dark grey G3

Yellow G1

Mustard G4

Orange A1

Red R1

Burgundy B2

Violet V4

Dark blue B3

Blue B1

Dark brown M1

Earth M2

Green V1

Forest green V2

Leaf green V3


Black N1


MARTE 2 | NEW PROFESSIONAL

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1121A30C1SB1 + Accessories Code 1218020







Model shown Code 1121A30R0SB1
+ Accessories Code 1218020



Model shown Code 1121A30C1CB1
+ Accessories Code 1218020

PARAMETRIC CODE

1	1	2	1	A	3				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm — 0 — 7

80 cm — 1 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

BREATHING HOLE

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

See colours on page 4

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 47-92 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +65° -50°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1218020
Collapsible cotsides	1221070
Safety magnet	1221030

ACCESSORIES	CODE
Safety double click	1221040
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

MARTE 3 | NEW PROFESSIONAL

CE

3-SECTION ELECTRICAL OR HYDRAULIC COUCH





Model shown Code 1121B30C1SB1 + Accessories Code 1218020






Sitting position







PARAMETRIC CODE

1	1	2	1	B	3				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm — 0 — 7

80 cm — 1 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PEDAL 2 ACTUATORS — P2 — 8

ALL AROUND 2 ACTUATORS — C2 — 8

BREATHING HOLE

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

See colours on page 4

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 47-92 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring+65° -50°
Middle section adjustment	with gas spring or electrically +20°
Footrest adjustment	with gas spring +85°
Trendelenburg position	+20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	80 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1218020
Safety magnet	1221030
Safety double click	1221040

ACCESSORIES	CODE
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060


ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3


MARTE 4 | NEW PROFESSIONAL

4-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1121C30C1SB1 + Accessories Code1218020







All around foot bar switch



Model shown Code 1121C30C1CB1
+ Accessories Code 1218020

PARAMETRIC CODE

1	1	2	1	C	3				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

80 cm

1

7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

ALL AROUND 1 ACTUATOR

C1

8

HYDRAULIC 0 ACTUATORS

R0

8

BREATHING HOLE

STANDARD

S

9

COMFORT

C

9

UPHOLSTERY

10

See colours on page 4

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 47-92 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +65° -50°
Armrests on headrest adjustment	manual
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	75 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1218020
Collapsible cotsides	1221070
Safety magnet	1221030

ACCESSORIES	CODE
Safety double click	1221040
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

MARTE 5 | NEW PROFESSIONAL

CE

5-SECTION ELECTRICAL OR HYDRAULIC COUCH





Model shown Code 1121D30C1SB1 + Accessories Code1218020






Trendelenburg position







PARAMETRIC CODE

1	1	2	1	D	3				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm — 0 — 7

80 cm — 1 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PEDAL 2 ACTUATORS — P2 — 8

ALL AROUND 2 ACTUATORS — C2 — 8

BREATHING HOLE

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

See colours on page 4

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 47-92 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +65° -50°
Armrests on headrest adjustment	manual
Middle section adjustment	with gas spring or electrically +20°
Footrest adjustment	with gas spring +85°
Trendelenburg position	+20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	85 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1218020
Safety magnet	1221030
Safety double click	1221040

ACCESSORIES	CODE
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060


ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

MONTELEONE | NEW PROFESSIONAL CE

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1121E30C1SB1 + Accessories Code 1218020





Hydraulic control



All around foot bar switch

PARAMETRIC CODE

1	1	2	1	E	3	0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

BREATHING HOLE

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

10

See colours on page 4

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 47-92 cm
Length	200 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +65° -50°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1218020
Safety magnet	1221030
Safety double click	1221040

ACCESSORIES	CODE
Swiss plug	1221050
Presa UK	1221060
UK 3-pin plug	1210060


ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3


SELENE 2 | NEW PROFESSIONAL

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1121F30C1CB1 + Accessories Code 1218020







Model shown Code 1121F30C1SB1 + Accessories Code 1218020



Retractable castors

PARAMETRIC CODE

1	1	2	1	F	3				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

80 cm

1

7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

ALL AROUND 1 ACTUATOR

C1

8

HYDRAULIC 0 ACTUATORS

R0

8

BREATHING HOLE

STANDARD

S

9

COMFORT

C

9

UPHOLSTERY

10

See colours on page 4

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 47-92 cm
Length	200 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +65° -50°
Trendelenburg Position	with gas spring +20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	75 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1218020
Safety magnet	1221030
Safety double click	1221040

ACCESSORIES	CODE
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060


ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

FEBO 6 | NEW PROFESSIONAL

CE

6-SECTION ELECTRICAL COUCH





Model shown Code 1121M30C3CG4 + Accessories Code 1218020





PARAMETRIC CODE

1	1	2	1	M	3	0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

ALL AROUND 3 ACTUATORS

C3

8

HAND CONTROL 3 ACTUATORS

H3

8

BREATHING HOLE

STANDARD

S

9

COMFORT

C

9

UPHOLSTERY

See colours on page 4

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical 58-102 cm
Length	210 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -15°
Armrests on headrest adjustment	manual
Middle section adjustment	electrically +70° -15°
Upper footrest adjustment	electrically +20° -15°
Lower footrest adjustment	with gas spring +40° -25°
Trendelenburg / AntiTrendelenburg position	with gas spring and electrically 15°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	105 kg

ACCESSORIES	CODE
Floor paper roll holder	1210020
Four castors	1218020
Safety magnet	1221030

ACCESSORIES	CODE
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

FEBO 6 PLUS | NEW PROFESSIONAL

6-SECTION ELECTRICAL COUCH




Model shown Code 1121N30C2SG4 + Accessories Code 1218020





All around foot bar switch



PARAMETRIC CODE

1	1	2	1	N	3	0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

ALL AROUND 2 ACTUATORS

C2

8

HAND CONTROL 2 ACTUATORS

H2

8

BREATHING HOLE

STANDARD

S

9

COMFORT

C

9

UPHOLSTERY

See colours on page 4

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical 58-102 cm
Length	210 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -15°
Armrests on headrest adjustment	manual
Middle section adjustment	electrically +70° -15°
Footrest adjustment	with gas spring +35° -15°
Trendelenburg / AntiTrendelenburg position	with gas spring and electrically 15°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	100 kg

OPTIONAL	CODICE
Floor paper roll holder	1210020
Four castors	1218020
Safety magnet	1221030

OPTIONAL	CODICE
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

FEBO 8 | NEW PROFESSIONAL

CE


8-SECTION ELECTRICAL COUCH





Model shown Code 1121P30C2SB1 + Accessories Code 1218020











PARAMETRIC CODE

1	1	2	1	P	3	0			
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

ALL AROUND 2 ACTUATORS

C2

8

HAND CONTROL 2 ACTUATORS

H2

8

BREATHING HOLE

STANDARD

S

9

COMFORT

C

9

UPHOLSTERY

10

See colours on page 4

TECHNICAL SPECIFICATIONS

Height adjustment	electrical 58-102 cm
Length	210 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -15°
Armrests adjustment	manual
Middle section adjustment	electrically +70° -15°
Footrests adjustment	with gas spring +35° -15°
Trendelenburg / AntiTrendelenburg position	with gas spring and electrically 15°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	24 cm
Under bed clearance (model with all round bar excluded)	12 cm
Weight	105 kg

ACCESSORIES	CODE
Floor paper roll holder	1210020
Four castors	1218020
Safety magnet	1221030

ACCESSORIES	CODE
Swiss plug	1221050
UK 3-pin plug	1221060
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 2 and 3

BREATHING HOLE

STANDARD

S





STANDARD

COMFORT

C





COMFORT

UPHOLSTERY COLOURS

Beige P1

Light grey G2

Taupe P2

Dark grey G3

Yellow G1

Mustard G4

Orange A1

Red R1

Burgundy B2

Violet V4

Dark blue B3

Blue B1

Dark brown M1

Earth M2

Green V1

Forest green V2

Leaf green V3

Black N1

FISIOTECH

SIMPLE

FISIOTECH



SIMPLE

EXAMINATION AND TREATMENT COUCHES

SIMPLE SERIES

Couches characterized by a frame with a basic design but especially thought to have clearance for patient hoists.

ELETTRA 1 | SIMPLE

2-SECTION ELECTRICAL OR HYDRAULIC COUCH





Model shown Code 1112A10R0SG3 + Accessories Code 1212010



Retractable castors







PARAMETRIC CODE

1	1	1	2	A	1				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001 → 1 (6)

WIDTH

65 cm → 0 (7)

80 cm → 1 (7)

MOVEMENT

PEDAL 1 ACTUATOR → P1 (8)

ALL AROUND 1 ACTUATOR → C1 (8)

HYDRAULIC 0 ACTUATORS → R0 (8)

PADDING

STANDARD → S (9)

COMFORT → C (9)

UPHOLSTERY

10 — See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 49-81 cm
Length	203 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring 0° +85°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	35 cm
Under bed clearance	14,5 cm
Weight	60 kg

ACCESSORIES	CODE
Four castors	1212010
Paper roll holder	1210010
Collapsible cotsides	1210030
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050

ACCESSORIES	CODE
Double-ply paper roll	1210060
Safety magnet	1212020
Safety double click	1212030
Swiss plug	1212040
UK 3-pin plug	1212041


ACCESSORIES see descriptions and pictures on page 106.


ELETTRA 2 | SIMPLE

2-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1112B10P1SB1 + Accessories Code 1212010







Comfort padding



Comfort padding

PARAMETRIC CODE

1	1	1	2	B	1				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1
6

WIDTH

65 cm
0
7

80 cm
1
7

MOVEMENT

PEDAL 1 ACTUATOR
P1
8

ALL AROUND 1 ACTUATOR
C1
8

HYDRAULIC 0 ACTUATORS
R0
8

PADDING

STANDARD
S
9

COMFORT
C
9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS

Height adjustment	electrical or hydraulic 49-81 cm
Length	190 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	35 cm
Under bed clearance	14,5 cm
Weight	60 kg

ACCESSORIES	CODE
Four castors	1212010
Paper roll holder	1210010
Collapsible cotsides	1210030
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050


ACCESSORIES	CODE
Double-ply paper roll	1210060
Safety magnet	1212020
Safety double click	1212030
Swiss plug	1212040
UK 3-pin plug	1212041


ACCESSORIES see descriptions and pictures on page 106.

ELETTRA 3 | SIMPLE

3-SECTION ELECTRICAL OR HYDRAULIC COUCH







Model shown Code 1112C10C1SG1 + Accessories Code 1212010



Sitting position



All around foot bar switch

PARAMETRIC CODE

1	1	1	2	C	1				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1
6

WIDTH

65 cm — 0 — 7

80 cm — 1 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PADDING

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-81 cm
Length	190 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +65° -65°
Middle section	fixed
Footrest adjustment	with gas spring +85°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	35 cm
Under bed clearance	14,5 cm
Weight	60 kg

ACCESSORIES	CODE
Four castors	1212010
Paper roll holder	1210010
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1212020
Safety double click	1212030
Swiss plug	1212040
UK 3-pin plug	1212041

ACCESSORIES see descriptions and pictures on page 106.

ELETTRA 4 | SIMPLE

4-SECTION ELECTRICAL OR HYDRAULIC COUCH




Model shown Code 1112D10P1SA1 + Accessories Code 1212010




Hydraulic control



Manual armrest adjustment





PARAMETRIC CODE

1	1	1	2	D	1				
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

65 cm — 0 — 7

80 cm — 1 — 7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PADDING

STANDARD — S — 9

COMFORT — C — 9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 49-81 cm
Length	190 cm
Width	65 / 80 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -20°
Armrests on headrest adjustment	manual
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	35 cm
Under bed clearance	14,5 cm
Weight	60 kg

ACCESSORIES	CODE
Four castors	1212010
Paper roll holder	1210010
Collapsible cotsides	1210030
Terry towel couch cover 65 cm	1210040
Terry towel couch cover 80 cm	1210050

ACCESSORIES	CODE
Double-ply paper roll	1210060
Safety magnet	1212020
Safety double click	1212030
Swiss plug	1212040
UK 3-pin plug	1212041

ACCESSORIES see descriptions and pictures on page 106.

FISIOTECH

BOBATH

FISIOTECH




BOBATH

COUCHES SUITABLE FOR
BOBATH CONCEPT

BOBATH SERIES


Large couches suitable for treatment with Bobath concept.




BOBATH 1 | BOBATH

1-SECTION ELECTRICAL COUCH







Model shown Code 1113A14P1NB1 + Accessories Code 1213010



Retractable castors



Collapsible cotsides

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 47-80 cm
Length	200 cm
Width	100 / 120 cm
Max. working load	200 kg
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	32 cm
Under bed clearance	17 cm
Weight	70 kg

ACCESSORIES	CODE
Four castors	1213010
Double foot pedal	1213030
Collapsible cotsides	1210030
Safety magnet	1213040

ACCESSORIES	CODE
Safety double click	1213050
Swiss plug	1213060
UK 3-pin plug	1213061

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

PARAMETRIC CODE

1	1	1	3	A	1		P1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

100 cm

4

7

120 cm

5

7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY


See colours and materials in the chart on the back cover flap.


10

BOBATH 2 PLUS | BOBATH

2-SECTION ELECTRICAL COUCH







Model shown Code 1113B15P1NV1 + Accessories Code 1213010



Headrest adjustment with gas spring



Retractable castors

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 47-80 cm
Length	200 cm
Width	100 / 120 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	32 cm
Under bed clearance	17 cm
Weight	70 kg

ACCESSORIES	CODE
Four castors	1213010
Double foot pedal	1213030
Collapsible cotsides	1210030
Safety magnet	1213040

ACCESSORIES	CODE
Safety double click	1213050
Swiss plug	1213060
UK 3-pin plug	1213061

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

PARAMETRIC CODE

1	1	1	3	B	1		P1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

100 cm

4

7

120 cm

5

7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY


See colours and materials in the chart on the back cover flap.


10

BOBATH 2 | BOBATH

2-SECTION ELECTRICAL COUCH







Model shown Code 1113C14P1NB1 + Accessories Code 1213010



Manual headrest adjustment



Retractable castors

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 47-80 cm
Length	200 cm
Width	100 / 120 cm
Max. working load	200 kg
Headrest adjustment	manual +60°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	32 cm
Under bed clearance	17 cm
Weight	70 kg

ACCESSORIES	CODE
Four castors	1213010
Double foot pedal	1213030
Collapsible cotsides	1210030
Safety magnet	1213040

ACCESSORIES	CODE
Safety double click	1213050
Swiss plug	1213060
UK 3-pin plug	1213061

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

PARAMETRIC CODE

1	1	1	3	C	1		P1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

100 cm — 4 — 7

120 cm — 5 — 7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

TITAN 1 | BOBATH

1-SECTION ELECTRICAL COUCH





Model shown Code 1113D15P1NB1 + Accessories Code 1213020





Retractable castors





PARAMETRIC CODE

1	1	1	3	D	1		P1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

100 cm — 4 — 7

120 cm — 5 — 7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical 49-90 cm
Length	200 cm
Width	100 / 120 cm
Max. working load	350 kg
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	30 cm
Under bed clearance	15 cm
Weight	100 kg

ACCESSORIES	CODE
Four castors	1213020
Double foot pedal	1213030
Collapsible cotsides	1210030
Safety magnet	1213040

ACCESSORIES	CODE
Safety double click	1213050
Swiss plug	1213060
UK 3-pin plug	1213061


ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.


TITAN 2 | BOBATH

2-SECTION ELECTRICAL COUCH



Model shown Code 1113E15P1NB1 + Accessories Code 1213020







Headrest adjustment with gas spring



Frame with max working load of 350 kg

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 49-90 cm
Length	200 cm
Width	100 / 120 cm
Max. working load	350 kg
Headrest adjustment	with gas spring +60°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	30 cm
Under bed clearance	15 cm
Weight	100 kg

ACCESSORIES	CODE
Four castors	1213020
Double foot pedal	1213030
Collapsible cotsides	1210030
Safety magnet	1213040

ACCESSORIES	CODE
Safety double click	1213050
Swiss plug	1213060
UK 3-pin plug	1213061

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

PARAMETRIC CODE

1	1	1	3	E	1		P1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

100 cm — 4 — 7

120 cm — 5 — 7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.


10

TITAN 1L | BOBATH

1-SECTION ELECTRICAL COUCH



Model shown Code 1113F16P1NB1 + Accessories Code 1213021







Retractable castors

PARAMETRIC CODE

1	1	1	3	F	1	6	P1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

180 cm

6

7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 49-90 cm
Length	200 cm
Width	180 cm
Max. working load	250 kg
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	30 cm
Under bed clearance	15 cm
Weight	100 kg

ACCESSORIES	CODE
Four castors	1213021
Double foot pedal	1213030
Collapsible cotsides	1210030
Safety magnet	1213040

ACCESSORIES	CODE
Safety double click	1213050
Swiss plug	1213060
UK 3-pin plug	1213061

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

ADE | BOBATH

1-SECTION FIXED-HEIGHT COUCH




Model shown Code 1113R15S0NV1




Adjustable foot



Model shown Code 1113R16S0NB1





TECHNICAL SPECIFICATIONS	
Fixed height	50 cm
Length	200 cm
Width	120 / 180 cm
Max. working load	200 kg
Adjustable foot	
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	50 kg

PARAMETRIC CODE

1	1	1	3	R	1		S0	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

120 cm

5

7

180 cm

6

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

ADE 2 | BOBATH

2-SECTION FIXED-HEIGHT COUCH




Model shown Code 1113S15S0NB1



Manual headrest adjustment





PARAMETRIC CODE

1	1	1	3	S	1	5	S0	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

120 cm

5

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

TECHNICAL SPECIFICATIONS	
Fixed height	50 cm
Length	200 cm
Width	120 cm
Max. working load	200 kg
Headrest adjustment	manual +60°
Adjustable foot	
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	50 kg

ACCESSORIES AND VARIANTS see descriptions and pictures on page 107.

FISIOTECH

TRACTION

FISIOTECH



TRACTION

TRACTION COUCHES

TRACTION SERIES

Specific couches for lumbar and cervical traction with or without traction unit.

ADONE | TRACTION

CE

FISIOTECH


3-SECTION ELECTRICAL COUCH



Model shown Code 1114A12P1SV1



Traction unit not supplied



SUPPLIED



TECHNICAL SPECIFICATIONS

Height adjustment	electrical 55-86 cm
Length	200 cm
Width	70 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -20°
Middle section	fixed
Main section	sliding by 12 cm
Chromed supports	as standard
With 4 castors	Ø 75 cm
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	32 cm
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Safety magnet	1214050
Safety double click	1214060

ACCESSORIES	CODE
Swiss plug	1214070
UK 3-pin plug	1214071

ACCESSORIES see descriptions and pictures on page 108.

PARAMETRIC CODE

1	1	1	4	A	1	2	P1	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

70 cm

2

7

MOVEMENT

PEDAL 1 ACTUATOR

P1

8

PADDING

STANDARD

S

9

UPHOLSTERY

10


See colours and materials in the chart on the back cover flap.


DAFNE | TRACTION


2-SECTION FIXED-HEIGHT COUCH



Model shown Code 1114B12S0NB1







SUPPLIED



PARAMETRIC CODE

1	1	1	4	B	1	2	S0	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

70 cm

2

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

TECHNICAL SPECIFICATIONS

Fixed height	65 cm
Length	200 cm
Width	70 cm
Max. working load	200 kg
Backrest tilt	10°
Main section	sliding by 12 cm
Adjustable feet	
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	70 kg

ACCESSORIES see descriptions and pictures on page 108.

SUPPLIED

ADONE - DAFNE

The accessories on this page are all included in Adone and Dafne couches.
You can reorder them separately, too, and specify when appropriate the upholstery colour.



Code 1415408
WEDGE CUSHION
dimensions 45x45x15 cm



Code 1415601
TRAPEZOID CUSHION
dimensions 40x50x30 cm



Code 1415201
HALF-ROLL CUSHION
dimensions 30x16x8 cm



Code 1214010
SMALL THORACIC HARNESS
min/max width 60/100 cm



Code 1214020
LARGE THORACIC HARNESS
min/max width 70/130 cm



Code 1214030
SMALL PELVIC HARNESS
min/max width 60/100 cm



Code 1214040
LARGE PELVIC HARNESS
min/max width 70/130 cm



Code 1216021
HEAD HALTER



Code 1216023
SPREAD BAR



Code 1214035
SPREADER BAR

FISIOTECH

TILT

FISIOTECH






TILT



TILT TABLES

TILT SERIES

Static supine therapy tables with electrically adjustable tilt up to 90°.



NETTUNO | TILT

1-SECTION TILT TABLE




Model shown Code 1115A30H1NB2





Axillary supports



PARAMETRIC CODE

1	1	1	5	A	3	0	H1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

HAND CONTROL 1 ACTUATOR

H1

8

PADDING

N

9

Without breathing hole.

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Fixed height	50 cm
Upholstery dimensions	190 x 65 cm
Overall dimensions	200 x 70 cm
Max. working load	150 kg
Angle adjustment	electrical 0° +90°
Complete with emergency lever for actuator unlocking	
Complete with four swivel castors, of which two locking	Ø 12,5 cm
Complete with safety straps	3
Footboard	fixed
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Under bed clearance	15 cm
Weight	70 kg

ACCESSORIES	CODE
Axillary supports	1215030
Padded patient restraint	1215040
10 cm safety strap	1215050
Patient table	1215060

ACCESSORIES	CODE
Safety magnet	1215070
Safety double click	1215080
Swiss plug	1215090
UK 3-pin plug	1215091

ACCESSORIES AND VARIANTS see descriptions and pictures on page 109.

GIASONE | TILT

1-SECTION TILT TABLE





Model shown Code 1115B30H1NB1





Independent footboards





TECHNICAL SPECIFICATIONS	
Fixed height	50 cm
Upholstery dimensions	190 x 65 cm
Overall dimensions	200 x 70 cm
Max. working load	150 kg
Angle adjustment	electrical 0° +90°
Complete with emergency lever for actuator unlocking	
Complete with four swivel castors, of which two locking	Ø 12,5 cm
Complete with safety straps	3
Independent footboards adjustment	manual +20° -20°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Under bed clearance	15 cm
Weight	70 kg

ACCESSORIES	CODE
Axillary supports	1215030
Padded patient restraint	1215040
10 cm safety strap	1215050
Patient table	1215060

ACCESSORIES	CODE
Safety magnet	1215070
Safety double click	1215080
Swiss plug	1215090
UK 3-pin plug	1215091

ACCESSORIES AND VARIANTS see descriptions and pictures on page 109.

PARAMETRIC CODE

1	1	1	5	B	3	0	H1	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

HAND CONTROL 1 ACTUATOR

H1

8

PADDING

N

9

Without breathing hole.

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

URANO | TILT

4-SECTION TILT TABLE





Model shown Code 1115C30H1SV1
+ Accessories Code 1215060





Patient table





PARAMETRIC CODE

1	1	1	5	C	3	0	H1	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

HAND CONTROL
1
ACTUATOR

H1

8

PADDING

STANDARD

S

9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Fixed height	50 cm
Upholstery dimensions	190 x 65 cm
Overall dimensions	200 x 70 cm
Max. working load	150 kg
Angle adjustment	electrical 0° +90°
Headrest adjustment	with gas spring +30° -80°
Complete with emergency lever for actuator unlocking	
Complete with four swivel castors, of which two locking	Ø 12,5 cm
Complete with safety straps	5
Complete with axillary supports	2
Independent footboards adjustment	manual +20° -20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Under bed clearance	15 cm
Weight	70 kg

ACCESSORIES	CODE
Padded patient restraint	1215040
10 cm safety strap	1215050
Patient table	1215060
Safety magnet	1215070

ACCESSORIES	CODE
Safety double click	1215080
Swiss plug	1215090
UK 3-pin plug	1215091

ACCESSORIES AND VARIANTS see descriptions and pictures on page 109.

EMERA | TILT

2-SECTION TILT TABLE





Model shown Code 1115D30H2SA1





Patient table





PARAMETRIC CODE

1	1	1	5	D	3	0	H2	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

HAND CONTROL
2 ACTUATORS

H2

8

PADDING

STANDARD

S

9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 54-96 cm
Upholstery dimensions	190 x 65 cm
Overall dimensions	200 x 70 cm
Max. working load	150 kg
Angle adjustment	electrical 0° +90°
Headrest adjustment	with gas spring +60° -20°
Complete with emergency lever for actuator unlocking	
Complete with four swivel castors, of which two locking	Ø 10 cm
Complete with safety straps	3
Complete with patient table	
Independent footboards adjustment	manual +20° -20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	15 cm
Weight	70 kg

ACCESSORIES	CODE
Axillary supports	1215030
Padded patient restraint	1215040
10 cm safety strap	1215050
Safety double click	1215070

ACCESSORIES	CODE
Safety double click	1215080
Swiss plug	1215090
UK 3-pin plug	1215091

ACCESSORIES AND VARIANTS see descriptions and pictures on page 109.

ULISSE | TILT

4-SECTION TILT TABLE









Model shown Code 1115F30H2SB1
+ Accessories Code 1215030 and 1215060

Patient table and axillary supports





TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 54-96 cm
Upholstery dimensions	190 x 65 cm
Overall dimensions	200 x 70 cm
Max. working load	150 kg
Angle adjustment	electrical 0° +90°
Headrest adjustment	with gas spring +30° -80°
Complete with emergency lever for actuator unlocking	
Complete with four swivel castors, of which two locking	Ø 10 cm
Complete with safety straps	5
Independent footboards adjustment	manual +20° -20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	none
Under bed clearance	15 cm
Weight	70 kg

ACCESSORIES	CODE
Axillary supports	1215030
Padded patient restraint	1215040
10 cm safety strap	1215050
Patient table	1215060

ACCESSORIES	CODE
Safety magnet	1215070
Safety double click	1215080
Swiss plug	1215090
UK 3-pin plug	1215091

ACCESSORIES AND VARIANTS see descriptions and pictures on page 109.

PARAMETRIC CODE

1	1	1	5	F	3	0	H2	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

65 cm

0

7

MOVEMENT

HAND CONTROL
2 ACTUATORS

H2

8

PADDING

STANDARD

S

9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

FISIOTECH


POSTURAL

FISIOTECH






POSTURAL



POSTURAL
REHABILITATION

POSTURAL SERIES


Couches designed for postural rehabilitation treatments.




POSEIDONE | POSTURAL CE

4-SECTION ELECTRICAL OR HYDRAULIC COUCH







Model shown Code 1117A39C1NB1 + Accessories Code 1218020





TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 56-100 cm
Length	200 cm
Width	56 cm
Max. working load	200 kg
Sliding lateral armrests rotating by	180°
Footrest adjustment	with gas spring +60° -90°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	100 (kg/m^3) / 1 cm
Lateral offset	none
Under bed clearance	12 cm
Weight	70 kg

ACCESSORIES	CODE
Height-adjustable pole + rope + ankle strap	1215010
Paper roll holder	1218010
Four castors	1218020
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1217010
Safety double click	1217020
Swiss plug	1217030
UK 3-pin plug	1217031

ACCESSORIES AND VARIANTS see descriptions and pictures on page 110.

PARAMETRIC CODE

1	1	1	7	A	3	9		N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

56 cm

9

7

MOVEMENT

PEDAL 1 ACTUATOR	ALL AROUND 1 ACTUATOR	ALUM. ALL AROUND 1 ACTUATOR	HYDRAULIC 0 ACTUATORS
P1	C1	A1	R0
8	8	8	8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

POSEIDONE 2 | POSTURAL CE


4-SECTION FIXED-HEIGHT COUCH




Model shown Code 1117E39S0NB1









TECHNICAL SPECIFICATIONS	
Fixed height	74 cm
Length	200 cm
Width	56 cm
Max. working load	180 kg
Sliding lateral armrests rotating by	180°
Footrest adjustment	with gas spring +60° -90°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	100 (kg/m^3) / 1 cm
Weight	60 kg

ACCESSORIES	CODE
Height-adjustable pole + rope + ankle strap	1215010
Paper roll holder	1218010

ACCESSORIES	CODE
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 110.

PARAMETRIC CODE

1	1	1	7	E	3	9	S0	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

56 cm

9

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY


See colours and materials in the chart on the back cover flap.


10

IRIS | POSTURAL

4-SECTION ELECTRICAL OR HYDRAULIC COUCH







Model shown Code 1117B10C1SB1 + Accessories Code 1211010





TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 45-90 cm
Length	190 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -20°
Sliding lateral armrests rotating by	180°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	17 cm
Weight	70 kg

ACCESSORIES	CODE
Height-adjustable pole + rope + ankle strap	1215010
Paper roll holder	1210010
Four castors	1211010
Double-ply paper roll	1210060

ACCESSORIES	CODE
Safety magnet	1217010
Safety double click	1217020
Swiss plug	1217030
UK 3-pin plug	1217031

ACCESSORIES AND VARIANTS see descriptions and pictures on page 110.

PARAMETRIC CODE

1	1	1	7	B	1	0		S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

65 cm

0

7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PADDING

STANDARD

S

9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

NIKE | POSTURAL


5-SECTION ELECTRICAL OR HYDRAULIC COUCH



Model shown Code 1117C10C1SN1 + Accessories Code 1211010







TECHNICAL SPECIFICATIONS	
Height adjustment	electrical or hydraulic 45-90 cm
Length	190 cm
Width	65 cm
Max. working load	200 kg
Headrest adjustment	with gas spring +60° -20°
Armrests adjustment	manual
Removable lateral armrests rotating by	180°
Middle section adjustment	with gas spring +20°
Footrest adjustment	with gas spring +75°
Trendelenburg position	+20°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Lateral offset	17 cm
Weight	70 kg

ACCESSORIES	CODE
Paper roll holder	1210010
Four castors	1211010
Double-ply paper roll	1210060
Safety magnet	1217010

ACCESSORIES	CODE
Safety double click	1217020
Swiss plug	1217030
UK 3-pin plug	1217031

ACCESSORIES AND VARIANTS see descriptions and pictures on page 110.

PARAMETRIC CODE

1	1	1	7	C	1	0		S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

65 cm

0

7

MOVEMENT

PEDAL 1 ACTUATOR — P1 — 8

ALL AROUND 1 ACTUATOR — C1 — 8

HYDRAULIC 0 ACTUATORS — R0 — 8

PADDING

STANDARD

S

9

UPHOLSTERY

10


See colours and materials in the chart on the back cover flap.

FISIOTECH

SET-CLOSE

FISIOTECH






SET-CLOSE



EXAMINATION
AND TREATMENT COUCHES

SET-CLOSE SERIES


Fixed-height and portable couches, with wooden or metal frame.



GEA | SET-CLOSE

2-SECTION FIXED-HEIGHT COUCH





Model shown Code 1116A10S0SG1



Paper roll holder



TECHNICAL SPECIFICATIONS	
Fixed height	78 cm
Length	190 cm
Width	65 cm
Max. working load	150 kg
Headrest adjustment	manual +60°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	50 kg

ACCESSORIES	CODE
Paper roll holder	1216010
Terry towel couch cover 65 cm	1210040

ACCESSORIES	CODE
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 112.

PARAMETRIC CODE

1	1	1	6	A	1	0	S0	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 9001

1

6

WIDTH

65 cm

0

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

STANDARD

S

9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

PONTO | SET-CLOSE


2-SECTION FIXED-HEIGHT COUCH




Model shown Code 1116B40S0SG2



Model shown Code 1116B43S0SV1





PARAMETRIC CODE

1	1	1	6	B	4		S0	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

CLEAR WOOD

4

6

WIDTH

65 cm

0

7

90 cm

3

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

STANDARD

S

9

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Fixed height	50 - 78 cm
Length	190 cm
Width	65 / 90 cm
Max. working load	130 kg
Headrest adjustment	manual +60°
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	50 kg

ACCESSORIES	CODE
Paper roll holder 65 cm	1216010
Paper roll holder 90 cm	1216020
Lower shelf 65 cm	1216040

ACCESSORIES	CODE
Lower shelf 90 cm	1216050
Terry towel couch cover 65 cm	1210040
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 112.

ARES | SET-CLOSE


2-SECTION FIXED-HEIGHT COUCH



Model shown Code 1116C40S0SB2







TECHNICAL SPECIFICATIONS	
Fixed height	78 cm
Length	190 cm
Width	65 cm
Max. working load	130 kg
Headrest adjustment	manual +60°
Tilting foldable legs	
Breathing hole and bung	as standard
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	50 kg

ACCESSORIES	CODE
Paper roll holder	1216010
Terry towel couch cover 65 cm	1210040

ACCESSORIES	CODE
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 112.

PARAMETRIC CODE

1	1	1	6	C	4	0	S0	S	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

CLEAR WOOD

4

6

WIDTH

65 cm

0

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

STANDARD

S

9

UPHOLSTERY


10

See colours and materials in the chart on the back cover flap.

REA | SET-CLOSE

1-SECTION FIXED-HEIGHT COUCH





Model shown Code 1116D43S0NB1



Lower shelf



TECHNICAL SPECIFICATIONS	
Fixed height	50 - 78 cm
Length	190 cm
Width	90 cm
Max. working load	130 kg
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	50 kg

ACCESSORIES	CODE
Paper roll holder	1216020
Lower shelf 90 cm	1216050

ACCESSORIES	CODE
Double-ply paper roll	1210060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 112.

PARAMETRIC CODE

1	1	1	6	D	4	3	S0	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

CLEAR WOOD

4

6

WIDTH

90 cm

3

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

N

9

Without breathing hole.

UPHOLSTERY

10

See colours and materials in the chart on the back cover flap.

HYDRA | SET-CLOSE

CE

FISIOTECH


WOODEN FRAME PORTABLE TABLE



Model shown Code 1116E40M0NB1







TECHNICAL SPECIFICATIONS	
Height	manual 60-75 cm
Length	218 cm
Width	65 cm
Dimensions closed	95 x 65 x 20 cm
Max. working load	150 kg
Headrest	removable
Armrest	removable
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 4 cm
Weight	17 kg

ACCESSORIES	CODE
HYDRA table carry case	1216070

ACCESSORIES AND VARIANTS see descriptions and pictures on page 112.

PARAMETRIC CODE

1	1	1	6	E	4	0	M0	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

CLEAR WOOD

4

6

WIDTH

65 cm

0

7

MOVEMENT

MANUAL COUCH

M0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY


10


See colours and materials in the chart on the back cover flap.

ANTERO | SET-CLOSE

2-SECTION STEEL FRAME PORTABLE TABLE







Model shown Code 1116F50S0NN1





PARAMETRIC CODE

1	1	1	6	F	5	0	S0	N	N1
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

CHROMED STEEL

5

6

WIDTH

64 cm

0

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

N1 Black


10

TECHNICAL SPECIFICATIONS	
Fixed height	72 cm
Length	187 cm
Width	64 cm
Dimensions closed	95 x 65 x 15 cm
Max. working load	95 kg
Headrest adjustment	manual +40°
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 2 cm
Weight	20 kg


ACCESSORIES	CODE
ANTERO table carry case	1216060

ACCESSORIES AND VARIANTS see descriptions and pictures on page 112.

FISIOTECH



GYNECOLOGY



FISIOTECH






GYNECOLOGY

EXAMINATION AND TREATMENT COUCHES

GYNECOLOGY SERIES

Couches for examination and treatment used in gynecology and urology.



ERA | GYNECOLOGY


3-SECTION ELECTRICAL GYNAECOLOGY ARMCHAIR




Model shown Code 1119A37H4NG1 + Accessories Code 1219030









TECHNICAL SPECIFICATIONS	
Height adjustment	electrical 53-85 cm
Length	180 cm
Width	54 cm
Max. working load	180 kg
Backrest adjustment	electrical 0° +90°
Leg section adjustment	electrical 0° -90°
Trendelenburg position	0° +21°
Armrests	folding
Leg supports adjustment	manual in height and rotation
Four locking castors	Ø 10 cm
Stainless steel bowl	removable
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 6 cm
Weight	72 kg

ACCESSORIES	CODE
Paper roll holder	1219020
Removable padded headrest	1219030
Double-ply paper roll	1210060
Safety magnet	1219040

ACCESSORIES	CODE
Safety double click	1219050
Swiss plug	1219060
UK 3-pin plug	1219061

ACCESSORIES AND VARIANTS see descriptions and pictures on page 113.

PARAMETRIC CODE

1	1	1	9	A	3	7	H4	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

54 cm

7

7

MOVEMENT

HAND CONTROL 4 ACTUATORS

H4

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY


See colours and materials in the chart on the back cover flap.


10

ESTIA | GYNECOLOGY

3-SECTION FIXED-HEIGHT GYNAECOLOGY COUCH







Model shown Code 1119B58S0NN1 + Accessories Code 1219010





Paper roll holder

PARAMETRIC CODE

1	1	1	9	B	5	8	S0	N	N1
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

CHROMED STEEL

5

6

WIDTH

62 cm

8

7

MOVEMENT

FIXED HEIGHT COUCH

S0

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

N1 Black

10

TECHNICAL SPECIFICATIONS	
Fixed height	76 cm
Length	185 cm
Width	62 cm
Max. working load	135 kg
Backrest adjustment	manual 0° +75°
Seat section	fixed
Footrest adjustment	manual 0° -90°
Leg supports adjustment	height and rotation
Stainless steel bowl	removable
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 5 cm
Weight	50 kg


ACCESSORIES	CODE
Paper roll holder	1219010

ACCESSORIES AND VARIANTS see descriptions and pictures on page 113.


EASY

FISIOTECH






EASY

BLOOD SAMPLING AND TREATMENT ARMCHAIR

EASY SERIES

Electrical blood sampling and treatment armchair, height-adjustable.



DEMETRA


BLOOD SAMPLING ARMCHAIR




Model shown Code 1120A37H4NB1 + Accessories Code 1219030









PARAMETRIC CODE

1	1	2	0	A	3	7	H4	N	
1	2	3	4	5	6	7	8	9	10

FRAME COLOUR

RAL 7035

3

6

WIDTH

54 cm

7

7

MOVEMENT

HAND CONTROL 4 ACTUATORS

H4

8

PADDING

Without breathing hole.

N

9

UPHOLSTERY

See colours and materials in the chart on the back cover flap.

10

TECHNICAL SPECIFICATIONS

Height adjustment	electrical 53-85 cm
Length	180 cm
Width	54 cm
Max. working load	180 kg
Backrest adjustment	electrical 0° +90°
Leg section adjustment	electrical 0° -90°
Trendelenburg position	0° +21°
Armrests	folding
Four locking castors	Ø 10 cm
Padding and upholstery fire-resistance rating	1IM
Padding: density/thickness	30 (kg/m^3) / 6 cm
Weight	72 kg

ACCESSORIES	CODE
Paper roll holder	1219020
Removable headrest	1219030
Double-ply paper roll	1210060
Safety magnet	1220010

ACCESSORIES	CODE
Safety double click	1220020
Swiss plug	1220030
UK 3-pin plug	1220031

ACCESSORIES AND VARIANTS see descriptions and pictures on page 114.

VARIANTI E OPTIONAL PER LETTINI

VARIATIONS AND OPTIONAL ACCESSORIES FOR COUCHES

VARIANTES ET OPTIONS POUR TABLES



FISIOTECH

FISIOTECH



UPHOLSTERY

TECHNICAL SPECIFICATIONS

FEATURE	TEST METHOD	UNIT	LIMIT VALUE
Tensile strength	UNI EN ISO 1421 Method 1	N	Lengthwise ≥ 210 Crosswise ≥ 90
Elongation strength	UNI EN ISO 1421 Method 1	%	Lengthwise ≥ 85 Crosswise ≥ 140
Determination of tear force	EN ISO 13937-2: 2000	daN	Lengthwise ≥ 2,5 Crosswise ≥ 2,2
Colour fastness to crocking	UNI EN ISO 105 X12: 2003 UNI EN ISO 105 D02: 1998	grey scale	Dry 4 Bleach 4/5 Ethyl alcohol 3/4
Determination of abrasion resistance: Martindale abrader	UNI EN ISO 5470-2: 2004	cycles	Method 1 (direct) 65.000 Method 2 (reverse) 210.000
Reaction ignition Fireproof class	UNI 9175: 1987 UNI 9175/FA1: 1994		CLASSE 1.IM
Free and hydrolyzed formaldehyde	UNI EN ISO 14184-1: 2011	mg/Kg	< 3
Phtalates	UNI EN ISO 14389: 2014	%	< 0,06
Arylamines	UNI EN ISO 14362-1: 2012	mg/Kg	Not detectable
Cleaning	Clean with damp soapy cloth and rinse well with clean water. Do not use solvents, bleaches, synthetic detergents, polishes or aerosol spray.		
Specific features	Antibacterial Fingicide		
Composition	Plasticised PVC 88% - Polyester 12%		


FIRE RESISTANCE ACCORDING TO:
IMO - 2010 FTP CODE Part 8
"Test for upholstered furniture"
par 7.2 SMOULDERING CIGARETTE TEST
par 7.3 PROPANE FLAME TEST
EN 1021-1: 2006 smouldering cigarette test
EN 1021-2: 2006 match flame equivalent test

UPHOLSTERY COLOURS

Beige P1


Light grey G2


Taupe P2


Dark grey G3


Yellow G1


Mustard G4


Orange A1


Red R1


Burgundy B2


Violet V4


Dark blue B3


Blue B1


Dark brown M1


Earth M2


Green V1


Forest green V2


Leaf green V3


Black N1

VARIATIONS

FOOT PEDAL
Standard foot pedal for height adjustment.



ALL AROUND FOOT BAR SWITCH
The height of the couch is adjusted by simply swinging the bar placed all around the base of the couch.



ALUMINIUM ALL AROUND FOOT BAR SWITCH
The height of the couch is adjusted by simply pressing the bar placed all around the base of the couch.



HYDRAULIC CONTROL
The height of the couch is adjusted by simply pressing the hydraulic piston pedal.



ACCESSORIES

PAPER ROLL HOLDER - Code 1218010
Metal structure to be placed under the final section of the couch.



CASTORS - Code 1218020
Four retractable castors to be used only to move the couch Ø 7,5 cm.



SAFETY MAGNET - Code 1218030
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1218040
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.

NOT AVAILABLE ON MOVEMENT WITH ALUMINIUM ALL-AROUND FOOT BAR SWITCH.



SWISS PLUG - Code 1218050



Code 1218050

UK 3-PIN PLUG - Code 1215051



Code 1218051

NEW PROFESSIONAL

VARIATIONS AND ACCESSORIES

VARIATIONS

ACCESSORIES

FOOT PEDAL
Standard foot pedal for height adjustment.



ALL AROUND FOOT BAR SWITCH
The height of the couch is adjusted by simply swinging the bar placed all around the base of the couch.



HYDRAULIC CONTROL
The height of the couch is adjusted by simply pressing the hydraulic piston pedal.



PAPER ROLL HOLDER - Code 1210010

Metal structure to be placed under the final section of the couch.



CASTORS - Code 1218020
Four retractable castors to be used only to move the couch Ø7,5 cm




COLLAPSIBLE COT SIDES - Code 1221070
Cot sides ensuring the patient's safety while lying on the couch. Easy to use, with a simple manual mechanism.



SAFETY MAGNET - Code 1221030
Removing the magnet all electrical regulation functions are prevented. This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1221040
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1221050

UK 3-PIN PLUG - Code 1221060


Code 1221050


Code 1221060

SIMPLE

VARIATIONS AND ACCESSORIES

VARIATIONS

ACCESSORIES

FOOT PEDAL
Standard foot pedal for height adjustment.



ALL AROUND FOOT BAR SWITCH
The height of the couch is adjusted by simply swinging the bar placed all around the base of the couch.



HYDRAULIC CONTROL
The height of the couch is adjusted by simply pressing the hydraulic piston pedal.



PAPER ROLL HOLDER - Code 1210010
Metal structure to be placed under the final section of the couch.



CASTORS - Code 1212010
Four retractable castors to be used only to move the couch Ø 7,5 cm.



COLLAPSIBLE COT SIDES - Code 1210030
Cot sides ensuring the patient's safety while lying on the couch.
Easy to use, with a simple manual mechanism.



SAFETY MAGNET - Code 1212020
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1212030
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1212040

UK 3-PIN PLUG - Code 1212041


Code 1212040


Code 1212041

BOBATH

VARIATIONS AND ACCESSORIES

VARIATIONS

ACCESSORIES

FOOT PEDAL
Standard foot pedal for height adjustment.



CASTORS BOBATH - Code 1213010
Four retractable castors to be used only to move the couch Ø 7,5 cm.
For BOBATH S couches



CASTORS TITAN
Four retractable castors to be used only to move the couch Ø 7,5 cm.
For TITAN couches Code 1213020
For TITAN 1L couches Code 1213021



COLLAPSIBLE COT SIDES - Code 1210030
Cot sides ensuring the patient's safety while lying on the couch.
Easy to use, with a simple manual mechanism.



DOUBLE FOOT PEDAL - Code 1213030
The height of the couch is adjusted from both sides of the couch thanks to this accessory.



SAFETY MAGNET - Code 1213040
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1213050
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1213060

UK 3-PIN PLUG - Code 1213061


Code 1213060


Code 1213061

SAFETY MAGNET - Code 1214050
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1214060
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1214070


Code 1214030

UK 3-PIN PLUG - Code 1214071


Code 1214031

VARIATIONS

ACCESSORIES

HAND CONTROL
Hand control for tilt adjustment and for height adjustment in height-adjustable tilt tables.



AXILLARY SUPPORTS - Code 1215030
The axillary supports can be adjusted in height, width and depth.



PATIENT TABLE - Code 1215060
The patient table can be used when the tilt table is in a supine static position.
The patient table is designed so that the patient can lay its forearms on it.
Its height and depth can be manually adjusted.



PADDED PATIENT RESTRAINT - Code 1215040
Padded band designed to provide higher and more comfortable support to the patient.



SAFETY STRAP HEIGHT 10 cm - Code 1215050
Safety strap for an easy support of the patient.
Height 10 cm.



SAFETY MAGNET - Code 1215070
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1215080
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1215090

UK 3-PIN PLUG - Code 1215091


Code 1215090


Code 1215091

POSTURAL

VARIATIONS AND ACCESSORIES

FOOT PEDAL
Standard foot pedal for height adjustment.



ALL AROUND FOOT BAR SWITCH
The height of the couch is adjusted by simply swinging the bar placed all around the base of the couch.





ALUMINIUM ALL AROUND FOOT BAR SWITCH
Available only on POSEIDONE couch.
The height of the couch is adjusted by simply pressing the bar placed all around the base of the couch.



HYDRAULIC CONTROL
The height of the couch is adjusted by simply pressing the hydraulic piston pedal.





POSTURAL

VARIATIONS AND ACCESSORIES

PAPER ROLL HOLDER - Code 1218010 For POSEIDONE couch.
PAPER ROLL HOLDER - Code 1210010 For IRIS and NIKE couch.
Metal structure to be placed under the final section of the couch.



CASTORS - Code 1218020 For POSEIDONE couch Ø 7,5 cm.



CASTORS - Code 1211010 For IRIS and NIKE couch Ø 6 cm.
Four retractable castors to be used only to move the couch.



ADJUSTABLE POLE - Code 1215010
Removable pole, adjustable in height and equipped with rope and ankle strap.
Not suitable for NIKE couch.



SAFETY MAGNET - Code 1217010
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1217020
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1217030

UK 3-PIN PLUG - Code 1217031


Code 1217030


Code 1217031

SET-CLOSE

ACCESSORIES

PAPER ROLL HOLDER - Code 1216010 For 65 cm wide fixed-height couches.



PAPER ROLL HOLDER - Code 1216020 For 90 cm wide fixed-height couches.
Metal structure to be placed under the final section of the couch.



LOWER SHELF - Code 1216040 For 65 cm wide fixed-height couches.
LOWER SHELF - Code 1216050 For 90 cm wide fixed-height couches.
Support surface for couches.



ANTERO CARRY CASE - Code 1216060
ANTERO couch carry case.



HYDRA CARRY CASE - Code 1216070
HYDRA couch carry case.



GYNECOLOGY

VARIATIONS AND ACCESSORIES

VARIATIONS

HAND CONTROL
Hand control for sections and height adjustment.



ACCESSORIES

PAPER ROLL HOLDER COUCH ESTIA - Code 1219010
Metal structure to be placed under the backrest of the couch.



PAPER ROLL HOLDER ERA ARMCHAIR - Code 1219020
Metal structure to be placed under the backrest of the couch.



HEADREST ERA ARMCHAIR - Cod. 1219030
Removable padded headrest to give the patient higher comfort.



SAFETY MAGNET - Code 1219040
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1219050
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1219060

UK 3-PIN PLUG - Code 1219061


Code 1219060


Code 1219061

EASY

VARIATIONS AND ACCESSORIES

VARIATIONS

HAND CONTROL
Hand control for sections and height adjustment.



ACCESSORIES

PAPER ROLL HOLDER DEMETRA ARMCHAIR - Code 1219020
Metal structure to be placed under the backrest of the couch.



HEADREST DEMETRA ARMCHAIR - Cod. 1219030
Removable padded headrest to give the patient higher comfort.



SAFETY MAGNET - Code 1220010
Removing the magnet all electrical regulation functions are prevented.
This avoids any accidental activation of the actuators.

SAFETY DOUBLE CLICK WITH MAGNET - Code 1220020
A software that allows the movement of the actuators only through a double input of the control. This avoids any accidental activation of the actuators.



SWISS PLUG - Code 1220030

UK 3-PIN PLUG - Code 1220031


Code 1220030


Code 1220031

FLOOR PAPER ROLL HOLDER - Code 1210020
Paper roll holder for examination, massage therapy, and treatment couches.
Easy to move, powder coated metal structure.
Dimensions: length cm 71, width cm 25, height cm 81.



COUCH COVER - Code 1210040 For 65 cm wide couches.
COUCH COVER - Code 1210050 For 80 cm wide couches.
Terry towel couch cover with elasticated hem. It is suitable for any treatment and examination couch.
Washable, grey fabric.



DOUBLE-PLY PAPER ROLL - Code 1210060
Pack of 6 rolls.
Dimensions: width cm 60, length mt 65.




RIABILITAZIONE
RÉÉDUCATION
REHABILITATION
FISIOTECH


FISIOTECH

PARALLEL BARS

FOLDABLE PARALLEL BARS



3 mt steel handrails Code 1311010





3 mt wooden handrails Code 1311020



Powder coated steel structure. Telescopic, height-adjustable, anatomic handrails available in two versions - wooden or powder coated steel. The knobs prevent the loosening and movement of the handrails when used by patients.

TECHNICAL SPECIFICATIONS	
Handrails height adjustment	69-100 cm
Distance between handrails	60 cm
Base width	70 cm
Handrails length	300 cm
Single handrail maximum load	135 kg

MODELS	CODE
3 mt foldable parallel bars with metal handrails	1311010
3 mt foldable parallel bars with wooden handrails	1311020



Handrails height adjustment system.

ACCESSORIES not available for foldable parallel bars.

PARALLEL BARS

PARALLEL BARS



6 mt steel handrails Code 1311140



2 mt steel handrails Code 1311120



3 mt steel handrails Code 1311110



4 mt steel handrails Code 1311130



3 mt parallel bars wooden handrails Code 1311100



Handrails height adjustment system.

Powder-coated steel structure with wooden platform with nonslip covering and convenient access ramps. Telescopic, height- and width-adjustable, anatomic handrails available in two versions - wooden or powder coated steel. The knobs prevent the loosening and movement of the handrail when used by patients.

TECHNICAL SPECIFICATIONS	
Handrails height adjustment	65-100 cm
Distance between handrails adjustment	40-80 cm
Platform width	59 cm
Base width	70 cm
Handrails length	200 / 300 / 400 / 600 cm
Single handrail maximum load	135 kg

MODELS	CODE
3 mt parallel bars with wooden handrails	1311100
3 mt parallel bars with steel handrails	1311110
2 mt parallel bars with steel handrails	1311120
4 mt parallel bars with steel handrails	1311130
6 mt parallel bars with steel handrails	1311140


ACCESSORIES see descriptions and pictures on page 122.

PARALLEL BARS

CE

FISIOTECH

CORNER PARALLEL BARS



Corner parallel bars with steel handrails Code 1311150

Powder-coated steel structure with wooden platform with nonslip covering and convenient access ramps. Telescopic height adjustable, anatomic powder coated steel handrails. The knobs prevent the loosening and movement of the handrail when used by patients.



Handrails height adjustment system.

TECHNICAL SPECIFICATIONS	
Handrails height adjustment	65-100 cm
Platform width	59 cm
Base width	70 cm
Length of the two platforms	200-270 cm
Single handrail maximum load	135 kg

MODELS	CODE
Corner parallel bars with steel handrails	1311150

ACCESSORIES see descriptions and pictures on page 122.

PARALLEL BARS

ACCESSORIES

PATIENT STABILISER - Code 1311200

The patient stabiliser may be used on all parallel bars with platform to support the patient during exercises. The patient stabiliser is equipped with four stoppers mountable on the ends of the handrails and a joining plate to connect a second patient stabilizer in case it is purchased. The double patient stabiliser helps the patient when reversing the direction.



Patient stabiliser wooden handrails



Double patient stabiliser wooden handrails





Please find below a series of varnished wood accessories that may be used to create gait-training courses within a specific therapy and rehabilitation programme.
They are only suitable on parallel bars with platform.



OBSTACLES / DIVIDER BOARD - Code 1311221

The obstacle / divider board set is made of eight varnished wood obstacles with magnet.
Dimensions 65 x 13 x 1,5 cm





CONCAVE / CONVEX PLATFORMS - Code 1311230

Made of two parts creating a concave or convex platform.
Dimensions 290 x 46 x 8,5 cm

Not suitable for parallel bars code 1311120 and 1311150.





PARALLEL BARS

ACCESSORIES

PLATFORM WITH DIVIDERS - Code 1311250
Create gait-training courses by putting accessories in the dividers.
Dimensions 290 x 50 x 3 cm

Not suitable for parallel bars code 1311120 and 1311150.



MOTOR TRAINING SANDBAGS - Code 1311260
Dimensions 38 x 20 cm



MOTOR TRAINING PVC BAGS - Code 1311270
Dimensions 38 x 20 cm



WOODEN ROLLER - Code 1311280
Length 20 cm Ø 3 cm



WOODEN STICK - Code 1311290
Length 20 cm Ø 1 cm



STAIRCASES

SINGLE FLIGHT REHABILITATION STAIRS



Single flight 4-step stairs Code 1312010



Inclined ramp optional Code 1312200

Powder coated steel structure, zinc-chromate sliding parts. Telescopic, height-adjustable, powder coated steel anatomic handrails. The knobs prevent the loosening and movement of the handrail when used by patients. Nonslip PVC covered closed steps.

TECHNICAL SPECIFICATIONS	
Handrails adjustment	manual 60-100 cm
Step width	56 cm
Step depth	30 cm
Step height – 4-step flight	13 cm
Step height – 3-step flight	16 cm
Single handrail maximum load	135 kg

MODELS	CODE
Single flight 4-step stairs	1312010
Single flight 3-step stairs	1312020

ACCESSORIES	CODE
Inclined ramp dimensions 170 x 49 cm	1312200



Handrails height adjustment system.

STAIRCASES

TWO FLIGHT REHABILITATION STAIRS



Stairs in corner position Code 1312030



Inclined ramp optional Code 1312200



Stairs in straight position Code 1312030

Powder coated steel structure, zinc-chromate sliding parts. Telescopic, height-adjustable, powder coated steel anatomic handrails. The knobs prevent the loosening and movement of the handrail when used by patients. Nonslip PVC covered closed steps in two different heights. The stairs can be set in a straight position as well as in a corner position.

TECHNICAL SPECIFICATIONS	
Handrails adjustment	manual 60-100 cm
Step width	56 cm
Step depth	30 cm
Step height – 4-step flight	13 cm
Step height – 3-step flight	16 cm
Single handrail maximum load	135 kg

MODELS	CODE
Two flight stairs	1312030

ACCESSORIES	CODE
Inclined ramp dimensions 170 x 49 cm	1312200



Handrails height adjustment system.

BALANCE TRAINING

ROCKER BOARDS

Varnished wood rocker boards with nonslip covering for balance training and rehabilitation.



Code 1313010

MODELS	CODE
With half-circle 170 x 60 x 14 cm (L x W x H) - Max working load 120 kg	1313010
With half-circle 80 x 60 x 13 cm (L x W x H) - Max working load 120 kg	1313020
With half-circle 60 x 40 x 10 cm (L x W x H) - Max working load 120 kg	1313030
With half-circle 40 x 35 x 9 cm (L x W x H) - Max working load 120 kg	1313040
Central fulcrum 60 x 40 x 10 cm (L x W x H) - Max working load 120 kg	1313050



Code 1313020 -1313030 -1313040



Code 1313050

BOBATH WOBBLE BOARDS

Varnished wood Bobath wobble boards with nonslip covering for balance training and rehabilitation.
With central semisphere



Code 1313060 -1313070

MODELS	CODE
Bobath wobble boards - Diameter 60 cm - Height 13 cm - Max load 120 kg	1313060
Bobath wobble boards - Diameter 80 cm - Height 13 cm - Max load 120 kg	1313070



SET OF THREE STEPS

Clear varnished wood steps with nonslip covering. They are suitable for all postural training therapies and are useful to promote proper musculoskeletal function, to work on body alignment, balance and coordination.
Two models with different heights - BIG and SMALL - are available.



Code 1313075 -1313076

MODELS	CODE
SMALL 60 x 60 x 10/15/20 cm (L x W x H) - Max load 180 kg	1313075
BIG 60 x 60 x 30/40/50 cm (L x W x H) - Max load 180 kg	1313076



BALANCE TRAINING

MOVIN' STEP - Code 1313151

It can either be used as a step for low-impact cardiovascular training or as a moving base to improve balance and co-ordination while exercising.
An integrated valve enables to adjust the air flow and to choose the correct movement range.



Code 1313151

TECHNICAL SPECIFICATIONS	
Length	35 cm
Width	35 cm
Height	7 cm
Recommended max load	120 kg
Tested max load	200 kg



STABILITY WHEEL - Code 1313181

This proprioceptive footboard is a valid aid to fitness and rehabilitation to carry out different kinds of exercises to improve balance, muscle tone and posture.



Code 1313181

TECHNICAL SPECIFICATIONS	
Diameter	40 cm
Height	9 cm
Recommended max load	120 kg
Tested max load	200 kg



CORE BALANCE - Code 1313171

It is made of a hard base and a soft green dome, which is removable from the base and inflatable according to individual needs. It is a multi-functional tool, which can be used for different kinds of training, from fitness to post-traumatic knee and ankle rehab.



Code 1313171

TECHNICAL SPECIFICATIONS	
Base diameter	47 cm
Base height	25 cm
Dome diameter	40 cm
Dome height	20 cm
Recommended max load	120 kg
Tested max load	200 kg



BALANCE TRAINING

DISCO SPORT - Code 1313191
Ideal for balance and motricity exercises, and for functional skill recovery.


Code 1313191

TECHNICAL SPECIFICATIONS	
Diameter	55 cm
Recommended max load	120 kg
Tested max load	200 kg
Colour	orange



DISCO DOME - Code 1313201
Ideal for step and balance exercises. The Disco Dome is made of a hard base and a soft orange disc, which can also be used separately for those who wish to perform more intense balance exercises.


Code 1313201

TECHNICAL SPECIFICATIONS	
Diameter	60 cm
Recommended max load	120 kg
Tested max load	200 kg



NEEDLE PUMP - Code 1411160
With flexible extension.


Code 1411160

MANUAL PUMP - Code 1411161
Ideal to inflate medium- and big-sized balls.


Code 1411161

ACCESSORIES







Accessory: Second tilting bar - Code 1314510

ARTEMIDE - Code 1314501

Rehabilitation armchair also suitable for upper and lower limb exercise. The barbell can be used on both sides. The weight/difficulty and angle are individually adjustable on various levels thanks to numerical markers. The legs can be blocked with straps.
Powder-coated steel structure. A second tilting bar is available as an optional.

TECHNICAL SPECIFICATIONS	
Length	150 cm
Width	60 cm
Height	70 -150 cm
Max working load	130 kg
Backrest adjustment	with gas spring 0° + 90°
Adjustable feet	4
Weights included	14 kg
Upholstery fire-resistance rating	1IM
Padding: density/thickness	40 (kg/m^3) / 4 cm
Weight	85 kg
Upholstery colour	blue

PEDAL EXERCISER FOR REHABILITATION - Code 1314520

Pedal exerciser to recover coordination and upper and lower limb muscle tone. Painted steel foldable structure, four non-slip tips, pedals with adjustable foot straps, handle for effort adjustment.

TECHNICAL SPECIFICATIONS	
Dimensions	49 x 37 x 34 cm
Minimum height pedals in horizontal position	27 cm
Maximum height pedals in vertical position	34 cm
Weight	2,3 kg



Code 1314520



With floor base Code 1315050



Wall-mounted Code 1315040

SHOULDER WHEEL

Shoulder and elbow rehabilitation instrument. The wheel is adjustable in height. Powder-coated steel structure with chrome-plated sliding part. Resistance mechanism for effort adjustment.

TECHNICAL SPECIFICATIONS	
Height	160 cm
Width	110 cm
Wall-mounted wheel depth	30 cm
Floor-mounted wheel depth	66 cm
Base of the wheel with base	90 x 60 cm

MODELS	CODE
Wall-mounted	1315040
With floor base	1315050

HAND REHABILITATION LADDER - Code 1315070

Wooden ladder for hand rehabilitation exercises.

TECHNICAL SPECIFICATIONS	
Width	30 cm
Height	60 cm



Code 1315070


SPIRALS – WORM SCREW

Vertical spiral for hand and arm torsion, pronation and supination.
Horizontal spiral for hand and forearm pronation and supination.
Worm screw for hand rotation.



Code 1315080

TECHNICAL SPECIFICATIONS	
Wooden base length	40 cm
Wooden base width	17 cm
Vertical spiral height	30 cm
Horizontal spiral height	25 cm
Worm screw height	20 cm



Code 1315090

MODELS	CODE
Vertical spiral	1315080
Horizontal spiral	1315090
Worm screw	1315100



Code 1315100

ERGOTHERAPY BOARD - Code 1315060

Upper limb rehabilitation board. Wall mountable, with nine elements simulating everyday activities like:

1. screwing/unscrewing a light bulb
2. turning on and off a switch
3. plugging/unplugging
4. sliding a bolt lock
5. rotating a latch
6. sliding a sliding lock
7. hanging a door chain
8. opening/closing a handle
9. opening/closing a handgrip



Code 1315060

TECHNICAL SPECIFICATIONS	
Length	100 cm
Width	60 cm

LIMB REHABILITATION

HANDGRIPS - Code 1315110
Instrument for hand and wrist rehabilitation. Sponge grip.
Pack of two.



Code 1315110

FOAM BALLS - Code 1315120
Soft, light and colorful foam ball, ideal for functional rehabilitation of fingers and hands.
Dimensions Ø 7 cm.




Code 1315120

WOODEN SPHERES - Code 1315130
Wooden spheres, ideal for functional rehabilitation of fingers and hands.
Set of four pieces. Assorted colours
Dimensions Ø 6 cm.



Code 1315130

WOODEN DISKS - Code 1315140
Colorful wooden disks, ideal for functional rehabilitation of fingers and hands.
Set of four pieces. Assorted colours
Dimensions Ø 6x3 cm.



Code 1315140

GAS PEDAL FOOTBOARD - Code 1315150
Made of varnished wood, with a leatherette soft wedge.
Suitable for ankle and foot exercises.



Code 1315150

TECHNICAL SPECIFICATIONS	
Length	24 cm
Width	14 cm
Height	15 cm





Fisiotek Code 1314010 - 1314020







FISIOTEK 3000

FISIOTEK 3000 series for passive rehabilitation of lower limbs. They are intended for use in rehabilitation to restore joint movement in both surgically and non-surgically treated medical conditions.

TECHNICAL SPECIFICATIONS	FISIOTEK 3000GS	FISIOTEK 3000TS
Hip and knee mobilization	•	•
Ankle mobilization	-	•
Speed control extension	•	•
Speed control flexion	•	•
Force control	•	•
Working time control	•	•
Extension automatic increase	•	•
Flexion automatic increase	•	•
Pause during extension	•	•
Pause during flexion	•	•
Warm Up cycles	• (*)	• (*)
Use memory card	•	-
Adjustable foot support	•	-
Extension automatic increase limit	•	-
Flexion automatic increase limit	•	-
Repetitions to extension limit	•	-
Repetitions to flexion limit	•	-
Knee range of movement	-10° ÷ 120° (*)	-10° ÷ 110°
Ankle range of movement	-7° ÷ 115°	-7° ÷ 115°
Hip range of movement	-	20° ÷ 0° ÷ 40°
Speed	0,8°÷ 3,5°/sec.	0,8°÷ 3,5°/sec.
Force	0 ÷ 40 kg	0 ÷ 40 kg
Weight	9,5 kg	14 kg

(*) Note: with the provided button.

MODELS	CODE
FISIOTEK 3000GS	1314010
FISIOTEK 3000TS	1314020

ACCESSORIES see descriptions and pictures on page 134.

LIMB REHABILITATION

ACCESSORIES FISIOTEK 3000GS - 3000TS

MOBILE PROGRAMMING KEYBOARD - Code 1314030
With its graphic, user-friendly display, the keyboard makes the equipment fully and simply programmable. The graphic symbols provided are self-explanatory.



Code 1314030




TROLLEY FOR INDOOR - Code 1314040
Easy and functional, it is fitted with nonslip supports for the Fisiotek machine and with a shelf.


Code 1314040

ACCESSORY FOR LIMBS SHORTER THAN 72 CM - Code 1314050
This accessory can only be mounted on FISIOTEK 3000 TS. It can be used in the rehabilitation of limbs between 61 cm and 72 cm in length for a range from 0° to 110°. Using same accessory, you can reach 135° of bending in a limb from 72 to 100 cm in length.


Code 1314050



Fisiotek LT-5R Code 1314060



FISIOTEK LT-5R - Code 1314060
Fisiotek LT-5R is used for the passive rehabilitation of the shoulder. It is used in rehabilitation to passively restore joint movement in both surgically and non-surgically treated medical conditions. Adjustable in height.



TECHNICAL SPECIFICATIONS	
MOVEMENTS CARRIED OUT	
Elevation in flexion range	0° ÷ 180°
Elevation in abduction	complete physiological range
Internal-external rotation range	90° ÷ 0° ÷ 90°
PROGRAMMABLE FUNCTIONS	
Speed	2° ÷ 4,5°/sec
Range of movement	fast electromechanical limit setting

ACCESSORIES

ARM SUPPORT - Code 1314070
Arm support for internal-external rotation movement.


Code 1314070

BATTERIES PACK AND BATTERY CHARGER - Code 1314080



Code 1314080

PULLEY THERAPY



Pulley therapy cage Code 1317001 (couch not included)

PULLEY THERAPY CAGE - Code 1317001
Self-supporting steel structure with 5 mm zinc-chromed steel wire mesh forming a 5x5 cm grid. Base with adjustable feet. Cage only.

TECHNICAL SPECIFICATIONS	
Max safe working load	300 kg
Weight of the structure	125 kg

ACCESSORIES

KIT FOR PULLEY THERAPY - Code 1317100
Kit with all the accessories for pulley therapy, slings, handles, weights and other useful tools.

No. 1 heel sling
No. 1 pair of ankle slings
No. 1 leg sling
No. 1 strap
No. 1 ankle sling
No. 1 leg sling strap
No. 1 arm sling strap
No. 1 head sling
No. 1 pelvis sling
No. 1 trunk sling
No. 1 heel sling
No. 1 20 mt long rope
No. 2 maniglie
No. 8 carrucole
No. 8 moschettoni
No. 2 cavigliere con gancio
No. 6 sandbags (one for each weight available)
No. 8 S hooks

Complete list and information on page page XX.

PULLEY THERAPY

ACCESSORIES

The accessories on this page are all included in the KIT FOR PULLEY THERAPY - Code 1317100.
You can also buy them separately.



Code 1317101
HEEL SLING



Code 1317102
PAIR OF ANKLE SLINGS



Code 1317103
LEG SLING



Code 1317104
STRAP



Code 1317105
ANKLE SLING



Code 1317106
LEG SLING STRAP



Code 1317107
ARM SLING STRAP



Code 1317108
HEAD SLING



Code 1317109
PELVIS SLING



Code 1317110
TRUNK SLING



Code 1317111
HEEL SLING WITH RINGS



Code 1317112
20 MT LONG ROPE

PULLEY THERAPY

ACCESSORIES

The accessories on this page are all included in the KIT FOR PULLEY THERAPY - Code 1317100.
You can also buy them separately.



Code 1317113
HANDLE



Code 1317114
PULLEY



Code 1317115
SNAP HOOK



Code 1316052
ANKLE STRAP WITH RING



Code 1317116
S HOOK



Code 1316031
0,5 kg sandbag with eyelet

Code 1316032
1 kg sandbag with eyelet

Code 1316033
2 kg sandbag with eyelet

Code 1316034
3 kg sandbag with eyelet


Code 1316035
4 kg sandbag with eyelet

Code 1316036
5 kg sandbag with eyelet

PULLEY THERAPY



Single exercise pulley Code 1315010



Double exercise pulley Code 1315020

SINGLE EXERCISE PULLEY - Code 1315010
Powder-coated steel single structure pulley for rehabilitation therapy. The working load can be adjusted up to 8 kg.

DOUBLE EXERCISE PULLEY - Code 1315020
Powder-coated steel structure with two pulleys with padded backrest upholstered with washable, fire-resistant leatherette. The working load can be adjusted up to 8 kg per side.

TECHNICAL SPECIFICATIONS	SINGLE	DOUBLE
Height	195 cm	185 cm
Width	30 cm	102 cm
Depth	28 cm	28 cm

ACCESSORIES

CAST IRON WEIGHT - Code 1316041
Varnished cast iron weight
With hole diameter 25 mm.



Code 1316041



Double exercise pulley with wall bar Code 1315030

DOUBLE EXERCISE PULLEY WITH WALL BAR - Code 1315030
Powder-coated steel structure with chrome-plated sliding parts. It is made of two lateral exercise pulleys with a central wall bar with varnished steel uprights and 13 wooden rungs. The working load can be adjusted up to 8 kg.

TECHNICAL SPECIFICATIONS	
Height	240 cm
Width	145 cm
Depth	28 cm

ACCESSORIES

CAST IRON WEIGHT - Code 1316041
Varnished cast iron weight
With hole diameter 25 mm.



Code 1316041


GINNASTICA
GYMNASTIQUE
GYMNASTICS
FISIOTECH


FISIOTECH

BALLS

GYMNIC CLASSIC

Ball used for various exercises in rehabilitation, in motor activity, for fitness and aerobics.
Used in the Pilates method and in functional, sports and athletic training, too.
Class I Medical Device.



TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	300 kg

MODELS	COLOURS	CODE
Ø 45 cm	Yellow	1411011
Ø 55 cm	Red	1411012
Ø 65 cm	Blue	1411013
Ø 75 cm	Yellow	1411014
Ø 85 cm	Red	1411015
Ø 95 cm	Blue	1411016
Ø 120 cm	Red	1411017



GYMNIC CLASSIC PLUS

The Gymnic Classic Plus is the same as Gymnic Classic, but with a greater safety for the user as in case of piercing, the ball deflates slowly. Class I Medical Device.



TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	300 kg



MODELS	COLOURS	CODE
Ø 55 cm	Red	1411021
Ø 65 cm	Blue	1411022
Ø 75 cm	Yellow	1411023



ACCESSORIES see descriptions and pictures on page 148.

BALLS

THERASENSORY

The Therasensory ball presents small bumps on its whole surface. This particular surface amplifies the sensory perceptions and provides, at the same time, a pleasant massage, by fostering blood circulation. Suitable for functional recovery and rehabilitation gymnastics.



TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	300 kg

MODELS	COLOURS	CODE
Ø 100 cm	Red	1411081
Ø 65 cm	Green	1411082



PHYSIO ROLL

Physio Roll was born from joining two balls together and creating an equipment with a great use range and a higher stability compared to balls. It used in the physiotherapy sector and it is ideal for those who have coordination and balance difficulties, for the disabled, for the elderly, and for nursery schools.



TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	300 kg



MODELS	COLOURS	CODE
Ø 30 cm - L 50 cm	Blue	1411071
Ø 40 cm - L 65 cm	Red	1411072
Ø 55 cm - L 90 cm	Yellow	1411073
Ø 70 cm - L 115 cm	Blue	1411074
Ø 85 cm - L 130 cm	Red	1411075



ACCESSORIES see descriptions and pictures on page 148.

BALLS

FOAM BALL - Cod. 1411031
Soft, light, and colourful foam ball ideal for muscular strength exercises, floor exercises and for lumbar support.
Assorted colours Dimensions: Ø 20 cm Weight: 0,150 kg


Code 1411031

SOFTGYM - Cod. 1411091
Soft, light, and resistant ball, which allows performing a wide range of exercises that vary from muscle strengthening to low-impact exercises. Foldable, so that it can be brought everywhere. Thanks to the straw included in the package, it is easily inflatable.
Colours: RED Dimensions: Ø 23 cm


Code 1411091

TECHNICAL SPECIFICATIONS	
Recommended max load	80 kg
Tested max load	100 kg



MEDICINE BALLS
Medicine balls are filled exclusively with water. Their surface allows achieving a good grip.



MODELS	COLOURS	CODE
Ø 10 cm - 500 gr	Green	1411051
Ø 12 cm - 1 kg	Red	1411052
Ø 15 cm - 2 kg	Yellow	1411053
Ø 17 cm - 3 kg	Blue	1411054
Ø 20 cm - 4 kg	Magenta	1411055
Ø 23 cm - 5 kg	Orange	1411056





ACCESSORIES see descriptions and pictures on page 148.

BALLS

RHYTHMIC BALLS

Balls suitable for training in rhythmic gymnastics.



MODELS	COLOURS	CODE
Max Ø 16 cm - 280 gr	Assorted colours	1411061
Max Ø 18,5 cm - 420 gr	Blue	1411062





SENSY-BALLS

Soft and pleasant to touch. Thanks to their bumps, they provide a healthy stimulation of the pressure points. The three available sizes allow a wide application range. The smaller version is particularly suitable for massage therapy and manipulative exercises, whereas the two larger versions are ideal for low-impact exercises. The three sizes are all suitable for reflexology and relaxing techniques.



TECHNICAL SPECIFICATIONS	
Recommended max load Ø 20 – Ø 28	80 kg
Tested max load	100 kg



MODELS		COLOURS	CODE
Ø 10 cm – 80 gr	Pack of two	Yellow	1411101
Ø 20 cm – 330 gr		Yellow	1411102
Ø 28 cm – 400 gr		Yellow	1411103



ACCESSORIES see descriptions and pictures on page 148.

REFLEX MASSAGE BALLS

The Reflexball reflex massage balls are ideal for massage, reflexology, hand therapy and relaxation exercises. Pack of two.



MODELS	COLOURS	CODE
Ø 8 cm - 85 gr	Green	1411111
Ø 9 cm - 120 gr	Blue	1411112
Ø 10 cm - 150 gr	Green	1411113




THERA FREEBALLS HAND - Code 1411131

Small balls that, thanks to their flexibility, allow performing specific exercises for hand rehabilitation.
Thera Freeballs Hand are ideal to develop finger mobility strengthening finger muscles, too. The set is made of two balls with different diameter.
Colours: RED and YELLOW Dimensions: 1 pc Ø 50 mm + 1 pc Ø 40 mm



Code 1411131



ACTIV ROLL - Code 1411121

A small roller with soft bumps for a comfortable massage. Activ Roll fosters blood circulation for a full-body benefit.
Used for foot and back massage, reflexology and skin-sensitisation exercises. The product is inflatable with a needle pump in order to achieve the desired level of inflation.
Colours: DARK GREY Dimensions: Ø 7,5 cm - L 15 cm



Code 1411121



ACCESSORIES see descriptions and pictures on page 148.

BALLS

ACCESSORIES

BALL SUPPORT "S" - Code 1411141
For balls Ø 45-55 cm.



BALL SUPPORT "M" - Code 1411151
For balls Ø 65-75 cm.



NEEDLE PUMP - Code 1411160
With flexible extension.



MANUAL PUMP - Code 1411161
Ideal for medium- or big-sized balls.



WALL BARS

WALL BARS WITH POWDER-COATED STEEL UPRIGHTS



Single section 13 rungs Code 1412010



Single section 14 rungs Code 1412020



Detail 14 rungs model

WALL BARS
Wall bars with oval clear-varnished beechwood rungs.
Powder-coated steel uprights with bars to fix the product to the wall.
One of the two upper rungs of the 14-rung wall bar is placed inwards while the other outwards to carry out suspension exercises.

TECHNICAL SPECIFICATIONS	
Height	240 cm
Width	90 cm
Maximum load	125 kg

MODELS	CODE
Single section 13 rungs	1412010
Single section 14 rungs	1412020

ACCESSORIES see descriptions and pictures on page 151.

WALL BARS

WOODEN WALL BARS



Single section Code 1412030



Double section Code 1412040

WOODEN WALL BARS
Clear-varnished beechwood wall bars with oval rungs.
One of the two upper rungs is placed inwards while the other outwards to carry out suspension exercises.
With zinc-plated steel joints to fix the wall bar to the wall.

TECHNICAL SPECIFICATIONS	
Height	240 cm
Maximum load	125 kg

MODELS	WIDTH	CODE
Single section	**90 cm**	1412030
Double section	**180 cm**	1412040

ACCESSORIES see descriptions and pictures on page 151.

WALL BARS

ACCESSORIES

CURVED BACK SUPPORT - Code 1412500

Curved back support for wall bars.
Padded, upholstered in black washable fire-resistant leatherette.
Colour: BLACK.

TECHNICAL SPECIFICATIONS	
Height	65 cm
Width	40 cm
Depth	18 cm



WALL BAR BOARD - Code 1412510

Wall bar board. Padded board upholstered in washable fire-resistant leatherette.
With feet-holding belt.
Colour: BLACK.

TECHNICAL SPECIFICATIONS	
Height	180 cm
Width	30 cm



SUSPENSION FRAME - Code 1412520

Suspension frame for wall bar. Zinc-chromed steel frame with grid for hooks.
Useful for cervical traction and pulley therapy.


TECHNICAL SPECIFICATIONS	
Height	40 cm
Width	64 cm
Depth	87 cm







MIRRORS



Screen-printed 100 x 200 mirror Code 1413010



Plain 100 x 200 mirror Code 1413020

SCREEN-PRINTED OR PLAIN MIRRORS

4mm-thick mirrors with accident prevention adhesive film. In case the mirror breaks, the fragments will not fall. Mirrors with safety film are mounted on a polystyrene board and a white MDF board. An aluminium frame binds together the components. The mirrors can either be mounted on a wall or on a moving inclination-adjustable trolley.

MODELS	CODE
Screen-printed 100 x 200 h cm. Grid 10 x 5 cm	1413010
Plain 100 x 200 h cm	1413020
Screen-printed 100 x 170 h cm. Grid 10 x 5 cm	1413030
Plain 100 x 170 h cm	1413040

ACCESSORIES

MIRROR TROLLEY

With powder-coated steel frame. When on the trolley, the mirror may be easily tilted and moved thanks to the four useful swivel castors placed on the base.

MODELS	CODE
Trolley for 100x200 cm mirrors	1413500
Trolley for 100x170 cm mirrors	1413510



MATTRESSES

FOAM MATTRESS

Blue PVC-covered compressed foam mattress. 600 gr/m2, antibacterial, anti-dust mite, fire-resistant class 2.
Removable upholstery, with zip and handles for easy transportation.
Colour: BLUE.

MODELS	DENSITY	CODE
200 x 100 x 5 cm	30 kg / m^3	1414011
200 x 100 x 10 cm	16 kg / m^3	1414012



AGGLOMERATED FOAM MATTRESS - Cod. 1414013

Agglomerated foam mattress with anti-mould antibacterial class 2 fire-resistant covering with non-slip bottom. Removable upholstery, with zip and handles for easy transportation.
Suitable for gym exercises.
Colour: BLUE.



Code 1414013

TECHNICAL SPECIFICATIONS	
Density	80 kg / m^3
Dimensions	200 x 100 x 5 cm



MAT GOLD® MATTRESS

Exclusive polyurethane integral skin foam formulation (without HFC and CFC) with ODP (Ozone Depletion Potential) blowing agents. The completely open cell structure gives the mat excellent shock absorption features, allowing the normal airflow from every single layer of the product, the homogeneity in thickness and density in any spot of the mat. The special skin on the top and on the bottom lends further mechanical resistance.
These mats guarantee great resistance to abrasions and an excellent grip.
Colour: BLUE.





MODELS	DENSITY	CODE
200 x 100 x 3 cm	PK14	1414031
200 x 100 x 4 cm	PK16	1414032
200 x 100 x 5 cm	PK20	1414033

HORIZONTAL MAT TROLLEY - Cod. 1414901

Varnished steel structure with four rubber swivel castors diameter 100 mm.
Dimensions 200x100 cm.



FOLDABLE MAT - Cod. 1414021
Foldable mat for floor exercises.
PVC upholstery.
Colour: BLUE Dimensions: 180 x 60 x 4 cm



Code 1414021

MAT - Cod. 1414022
Mat that can be easily rolled up.
Colour: GREY Dimensions: 180 x 50 x 1 cm



Code 1414022



FOLDABLE MAT - Cod. 1414023
Foldable mat for floor exercises.
Colour: GREY Dimensions: 180 x 50 x 1 cm



Code 1414023

AIREX®

GYM MATS



AIREX® GYM MATS



Airex AG is the world-leading manufacturer of gym mats under the brand AIREX®. Tested materials and proven Swiss quality in the manufacturing process guarantee product durability even under tough conditions. AIREX® products meet or exceed current legal regulations in regards to materials and accident prevention measures.

Reasons why you should choose AIREX® professional products

- Comfort to touch
- Simple to clean, antimicrobial finishing
- Hard-wearing material for many years of use
- Yielding and supportive, protection against injuries
- Excellent damping capacity and great stability in balance exercises
- Lightness and flexibility, easy to roll and store
- Water repellent. Closed cell foam. No penetration of water or dirt

Three tips for your AIREX® products

1. Use the products barefoot or in socks (no sneakers)
2. If your products get dirty, clean them with ordinary soapy water. Disinfection is not going to be necessary thanks to integrated Sanitized protection!
3. Let the AIREX® products air dry thoroughly if they get wet (they can also be used outdoors)

CORONELLA 185

AIREX® CORONELLA is suitable for: pilates, fitness, yoga, physiotherapy and rehabilitation, and club sport.

MODELS	COLOURS	WEIGHT	CODE
185 x 60 x 1,5 cm	Blue	2,2 kg	1414141
185 x 60 x 1,5 cm	Red	2,2 kg	1414142
185 x 60 x 1,5 cm	Slate	2,2 kg	1414143





Code 1414141



Code 1414142



Code 1414143

CORONELLA 200

AIREX® CORONELLA is suitable for: pilates, fitness, yoga, physiotherapy and rehabilitation, and club sport.

MODELS	COLOURS	WEIGHT	CODE
200 x 60 x 1,5 cm	Slate	2,4 kg	1414144
200 x 60 x 1,5 cm	Blue	2,4 kg	1414145
200 x 60 x 1,5 cm	Red	2,4 kg	1414146





Code 1414144



Code 1414145



Code 1414146

ACCESSORIES see descriptions and pictures on page 161.

AIREX®

GYM MATS

CORONA

AIREX® CORONA is suitable for: pilates, fitness, yoga, physiotherapy and rehabilitation and club sport.

MODELS	COLOURS	WEIGHT	CODE
185 x 100 x 1,5 cm	Blue	3,8 kg	1414151
185 x 100 x 1,5 cm	Red	3,8 kg	1414152
185 x 100 x 1,5 cm	Slate	3,8 kg	1414153





Code 1414151



Code 1414152



Code 1414153

ATLAS

AIREX® ATLAS is suitable for: fitness, physiotherapy and rehabilitation and club sport.

MODELS	COLOURS	WEIGHT	CODE
200 x 125 x 1,5 cm	Red	5,5 kg	1414161
200 x 125 x 1,5 cm	Waterblue	5,5 kg	1414162



Code 1414144



Code 1414145



ACCESSORIES see descriptions and pictures on page 161.

AIREX®

GYM MATS

HERCULES

AIREX® HERCULES is suitable for: physiotherapy and rehabilitation and club sport.

MODELS	COLOURS	WEIGHT	CODE
200 x 100 x 2,5 cm	Blue	6,6 kg	1414171
200 x 100 x 2,5 cm	Waterblue	6,6 kg	1414172





Code 1414171



Code 1414172

FITLINE 180

AIREX® FITLINE 180 is suitable for: fitness, pilates, yoga, and club sport.

MODELS	COLOURS	WEIGHT	CODE
200 x 58 x 1 cm	Waterblue	1,5 kg	1414181
200 x 58 x 1 cm	Charcoal	1,5 kg	1414182
200 x 58 x 1 cm	Kiwi	1,5 kg	1414183
200 x 58 x 1 cm	Pink	1,5 kg	1414184





Code 1414181



Code 1414182



Code 1414183




Code 1414184

ACCESSORIES see descriptions and pictures on page 161.

CORONELLA 120 - Code 1414191

AIREX® CORONELLA 120 is suitable for: fitness and club sport.


MODELS	COLOURS	WEIGHT	CODE
120 x 60 x 1,5 cm	Blue	1,4 kg	1414191



FITLINE 140 - Code 1414201

AIREX® FITLINE 140 is suitable for: fitness, pilates, yoga, and club sport.

MODELS	COLOURS	WEIGHT	CODE
140 x 58 x 1 cm	Waterblue	1,2 kg	1414201



PILATES 190 - Code 1414212

AIREX® PILATES 190 is suitable for: fitness, pilates, yoga.

MODELS	COLOURS	WEIGHT	CODE
190 x 60 x 0,8 cm	Anthracite	1,4 kg	1414212



ACCESSORIES see descriptions and pictures on page 161.

EYELETS

MODELS	CODE
Type 15 eyelets suitable for: AIREX® CORONELLA	1414301
Type 10 eyelets suitable for: AIREX® FITLINE140/180	1414302
Type 8 eyelets suitable for: AIREX® PILATES 190	1414303



WALL SUPPORT - Code 1414311
Wall support for AIREX® mats with eyelet.
Suitable for: CORONELLA – FITLINE140/180 – PILATES 190



WALL SUPPORT - Code 1414312
Wall support for AIREX® mats without eyelet.
Suitable for: CORONELLA – FITLINE 140/180



WALL SUPPORT - Code 1414313
Wall support for AIREX® mats without eyelet.
Suitable for: CORONA – HERCULES



AIREX®

BALANCE TRAINING

BALANCE PAD - Code 1414221

Smooth surface ideal for balance training. The thickness of 6 cm makes the product unstable. This, added to the yielding foam the pad is made of, causes a constant challenge to the body in maintaining balance and stabilising the joints.
Suitable for: pilates, fitness, yoga, physiotherapy and rehabilitation, club sport and balance training.



Code 1414221

MODELS	COLOURS	WEIGHT	CODE
41 x 50 x 6 cm	Blue	0,68 kg	1414221

BALANCE PAD ELITE

The AIREX® BALANCE PAD ELITE differs from the BALANCE PAD in that it has a honeycomb structure on the top and on the bottom.



Code 1414231



Code 1414232

MODELS	COLOURS	WEIGHT	CODE
41 x 50 x 6 cm	Blue	0,68 kg	1414231
41 x 50 x 6 cm	Kiwi	0,68 kg	1414232
41 x 50 x 6 cm	Lava	0,68 kg	1414233
41 x 50 x 6 cm	Pink	0,68 kg	1414234



Code 1414233



Code 1414234

BALANCE PAD XL - Code 1414251

The AIREX® BALANCE PAD XL is same as BALANCE PAD ELITE but with different dimensions.



Code 1414251

MODELS	COLOURS	WEIGHT	CODE
98 x 41 x 6 cm	Blue	1,4 kg	1414251

BALANCE BEAM - Code 1414241

Product developed especially for clinical balance and proprioception exercises. It is based on the shape of a conventional balance beam and used mainly in physiotherapy to regain walking and foot functions.



Code 1414241


MODELS	COLOURS	WEIGHT	CODE
160 x 24 x 6 cm	Blue	0,99 kg	1414241

CUSHIONS

CUSHIONS

Cushions available in different shapes and sizes.
Polyurethane foam padding density: 30 kg/m^3, Italian fire resistance rating: 1 IM.
Washable leatherette upholstery, Italian fire resistance rating: 1 IM.






ROLL

MODELS	COLOURS	CODE
L 20 x Ø 8 cm	assorted colours	1415101
L 35 x Ø 15 cm	assorted colours	1415102
L 50 x Ø 15 cm	assorted colours	1415103
L 30 x Ø 30 cm	assorted colours	1415104
L 50 x Ø 25 cm	assorted colours	1415105
L 100 x Ø 25 cm	assorted colours	1415106
L 100 x Ø 40 cm	assorted colours	1415107
L 60 x Ø 30 cm	assorted colours	1415108
L 120 x Ø 30 cm	assorted colours	1415109

See colours and materials in the chart on the back cover flap.





HALF-ROLL


MODELS	COLOURS	CODE
L 30 x Ø 16 x h 8 cm	assorted colours	1415201
L 50 x Ø 25 x h 12,5 cm	assorted colours	1415202
L 60 x Ø 30 x h 15 cm	assorted colours	1415203


See colours and materials in the chart on the back cover flap.

CUSHIONS

CUSHIONS

Cushions available in different shapes and sizes.
Polyurethane foam padding density: 30 kg/m^3, Italian fire resistance rating: 1 IM.
Washable leatherette upholstery, Italian fire resistance rating: 1 IM.






CUBE

MODELS	COLOURS	CODE
L 40 x P 40 x h 40 cm	assorted colours	1415301
L 30 x P 30 x h 30 cm	assorted colours	1415302

See colours and materials in the chart on the back cover flap.





PARALLELEPIPED

MODELS	COLOURS	CODE
L 40 x P 20 x h 10 cm	assorted colours	1415501
L 25 x P 50 x h 25 cm	assorted colours	1415502


See colours and materials in the chart on the back cover flap.

CUSHIONS

CUSHIONS

Cushions available in different shapes and sizes.
Polyurethane foam padding density: 30 kg/m^3, Italian fire resistance rating: 1 IM.
Washable leatherette upholstery, Italian fire resistance rating: 1 IM.






WEDGE

MODELS	COLOURS	CODE
L 35 x P 35 x h 10 cm	assorted colours	1415402
L 25 x P 25 x h 10 cm	assorted colours	1415403
L 45 x P 60 x h 15 cm	assorted colours	1415404
L 60 x P 60 x h 15 cm	assorted colours	1415405
L 60 x P 60 x h 30 cm	assorted colours	1415406
L 45 x P 45 x h 15 cm	assorted colours	1415408

See colours and materials in the chart on the back cover flap.





CUT-OUT WEDGE


MODELS	COLOURS	CODE
L 40 x P 34 x h 9 cm	assorted colours	1415407


See colours and materials in the chart on the back cover flap.

CUSHIONS

CUSHIONS

Cushions available in different shapes and sizes.
Polyurethane foam padding density: 30 kg/m^3, Italian fire resistance rating: 1 IM.
Washable leatherette upholstery, Italian fire resistance rating: 1 IM.





TRAPEZOID

MODELS	COLOURS	CODE
L 50 x P 40 x h 30 cm	assorted colours	1415601
L 50 x P 40 x h 30-16 cm	assorted colours	1415602

See colours and materials in the chart on the back cover flap.





MOVIN'SIT - Code 1417141

Movin' Sit is an inflatable dynamic cushion that activates the trunk muscles and supports the backbone. Its shape was carefully studied in order to promote good posture by facilitating the search for a balanced pelvic position.

TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	200 kg

MODELS	COLOURS	CODE
L 36 x P 36 x h 7 cm	Blue	1417141





DISC'O'SIT - Code 1417131
Disc'o'Sit sitting cushion for a correct posture. When positioned on a chair, it relieves the backbone, thanks to the movement it provides.
The Disc'o'Sit is useful in functional training as well.
Suitable for floor exercises to improve coordination, reflexes and stability.

TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	200 kg

MODELS	COLOURS	CODE
Ø 39 cm	Blue	1417131


ACCESSORIES



ACCESSORIES - Code 1417132
Cover with zip useful when the Disc'o'Sit is used as a sitting cushion.
Colours: BLUE

HOOPS


Colourful nylon hoops, flat section. Used for psychomotor courses along with cones with holes (code 1416031) and sticks (code 1416061).



MODELS	COLOURS	CODE
Ø 60 cm	assorted colours	1416011
Ø 70 cm	assorted colours	1416012
Ø 80 cm	assorted colours	1416013

PVC STICKS SET - Code 1416061

Set of four colourful gymnastics PVC sticks. Used for psychomotor courses along with cones with holes (code 1416031) and hoops (codes 1416011-12-13).



Code 1416061

MODELS	COLOURS	CODE
Ø 3 cm	assorted colours	1416061

CONE WITH HOLES - Code 1416031

Cone with holes available in four random colours (red, yellow, green, blue).
Split vertex to insert PVC sticks (code 1416061) and hoops (codes 1416011-12-13) and three symmetrical lateral holes to insert PVC sticks in order to create courses.



Code 1416031

MODELS	COLOURS	CODE
30 cm H	assorted colours	1416031

CLUBS AND STICKS

WOODEN CLUBS
Clear-varnished beech clubs.



MODELS	DIMENSIONS	WEIGHT	CODE
Small clubs	310 x 56 mm	300 gr	1416041
Medium clubs	350 x 56 mm	400 gr	1416042
Big clubs	390 x 56 mm	500 gr	1416043

NYLON CLUBS - Code 1416045
Pair of colourful nylon clubs for rhythmic gymnastics.



Code 1416045

MODELS	DIMENSIONS	WEIGHT	CODE
Nylon clubs	450 mm	150 gr	1416045

WOODEN STICKS
Clear-varnished wood stick for gymnastics.



LENGTH	DIAMETER	CODE
70 cm	Ø 2,4 cm	1416051
80 cm	Ø 2,4 cm	1416052
90 cm	Ø 2,4 cm	1416053
100 cm	Ø 2,4 cm	1416054
120 cm	Ø 2,4 cm	1416055
150 cm	Ø 2,4 cm	1416056



ACCESSORIES

WALL-MOUNTED STICK RACK - Code 1416071
Wall-mounted stick rack. Varnished steel tube.
Capacity: 100 pcs.
Dimensions: length 50 x width 35 x height 35 cm.



Code 1416071


ELASTICS

ELASTIC BANDS

Elastic bands are simple and practical tools for fitness, pilates, muscle toning and strengthening, rehabilitation and stretching.
Dimensions: 120 x 14,5 cm.



MODELS	COLOURS	CODE
Extra light	Yellow	1417011
Light	Red	1417012
Medium	Green	1417013
Heavy	Blue	1417014



25MT-LONG ELASTIC BANDS

Dimensions: 25 mt x 14,5 cm.



MODELS	COLOURS	CODE
Extra light	Yellow	1417021
Light	Red	1417022
Medium	Green	1417023
Heavy	Blue	1417024

ELASTIC BANDS WITH HANDLES

Polyurethane latex free elastic band rolls. Thanks to their dimensions, they are suitable for low-impact exercises and for use in pools, as they do tolerate chlorine.
Dimensions: 160 x 5 cm.



MODELS	COLOURS	CODE
Extra light	Red	1417031
Light	Yellow	1417032
Medium	Blue	1417033
Heavy	Green	1417034



ACCESSORIES

ELASTIC BAND HANDLES - Code 1417111



Code 1417111

ELASTICS

ROLL 20 mt x 5 cm

Polyurethane latex free elastic band rolls. Thanks to their dimensions, they are suitable for low-impact exercises and for use in pools, as they do tolerate chlorine.
Dimensions 20 mt x 5 cm.



MODELS	COLOURS	CODE
Extra light	Red	1417041
Light	Yellow	1417042
Medium	Blue	1417043
Heavy	Green	1417044

100% LATEX FREE

POWER LOOP

Power loop is a large elastic loop with a high resistance. It is ideal for functional training, fitness, physical therapy.
Thickness: 0,45 cm.



MODELS	DIMENSIONS	COLOURS	CODE
Extra light	104 x 1,25 cm	Red	1417051
Light	104 x 1,90 cm	Black	1417052
Medium	104 x 2,85 cm	Yellow	1417053
Strong	104 x 4,45 cm	Green	1417054
Extra strong	104 x 6,35 cm	Blue	1417055
Super extra strong	104 x 8,25 cm	Orange	1417056



TUBING 10 MT

Polyurethane, latex free high resistance tubing bands. They are a simple and practical tool for muscle toning and strengthening, for rehabilitation and stretching.
Dimensions: 10 mt.



MODELS	COLOURS	CODE
Extra light	Red	1417071
Light	Yellow	1417072
Medium	Violet	1417073
Heavy	Black	1417074

100% LATEX FREE

TUBING HANDLES - Code 1417112



Code 1417112

ELASTICS

TUBING BANDS WITH HANDLES

Polyurethane, latex free high resistance tubing bands. They are a simple and practical tool for muscle toning and strengthening, for rehabilitation and stretching.
Dimensions: 160 cm.

MODELS	COLOURS	CODE
Extra light	Red	1417061
Light	Yellow	1417062
Medium	Violet	1417063
Heavy	Black	1417064



EIGHT-SHAPED TUBING BANDS

Eight-shaped high resistance tubing bands with handles. Polyurethane, latex free elastics, ideal for functional training, fitness, physical therapy.

MODELS	COLOURS	CODE
Light	Red	1417101
Medium	Yellow	1417102
Strong	Violet	1417103
Extra strong	Green	1417104





TUBING BAND RINGS

Polyurethane, latex free high resistance tubing band rings. They are a simple and practical tool for muscle toning and strengthening, for rehabilitation and stretching.

MODELS	COLOURS	CODE
Light	Red	1417081
Medium	Yellow	1417082
Strong	Violet	1417083
Extra strong	Black	1417084



TUBING BAND RINGS WITH HANDLES

Polyurethane, latex free high resistance tubing band rings with handles ideal for functional training, fitness, physical therapy.

MODELS	COLOURS	CODE
Light	Red	1417091
Medium	Yellow	1417092
Strong	Violet	1417093
Extra strong	Black	1417094



ANKLE AND WRIST WEIGHTS

LEATHERETTE ANKLE AND WRIST WEIGHTS

Leatherette covering. Hook and loop fasteners.

MODELS	CODE
Pair 1 kg (2 x 0.5 kg each)	1417301
Pair 2 kg (2 x 1 kg each)	1417302
Pair 3 kg (2 x 1.5 kg each)	1417303
Pair 4 kg (2 x 2 kg each)	1417304
Pair 5 kg (2 x 2.5 kg each)	1417305
Pair 6 kg (2 x 3 kg each)	1417306
Pair 7 kg (2 x 3.5 kg each)	1417307
Pair 8 kg (2 x 4 kg each)	1417308



NEOPRENE ANKLE AND WRIST WEIGHTS

Neoprene covering. Hook and loop fasteners.
Colour: BLACK.

MODELS	CODE
Pair 1 kg (2 x 0.5 kg each)	1417311
Pair 2 kg (2 x 1 kg each)	1417312
Pair 3 kg (2 x 1.5 kg each)	1417313
Pair 4 kg (2 x 2 kg each)	1417314
Pair 5 kg (2 x 2.5 kg each)	1417315







DUMBBELLS



DUMBBELLS
Cast iron dumbbells covered in vinyl for a good grip.
Pack of two.

MODELS	CODE
1 kg dumbbells	1417401
2 kg dumbbells	1417402
3 kg dumbbells	1417403
4 kg dumbbells	1417404
5 kg dumbbells	1417405

ACCESSORIES

DUMBBELL RACK - Code 1417451



Code 1417451

FIT-BALL ROLLER - Code 1418121

The Fit-Ball Roller is an effective alternative to hard cylinders. This multi-purpose tool allows performing a wide range of exercises, which vary from abdominal strengthening to stretching, relaxation and Pilates. Inflatable according to the desired hardness.



Code 1418121

TECHNICAL SPECIFICATIONS	
Recommended max load	120 kg
Tested max load	200 kg
Diameter	18 cm
Length	75 cm



FOAM-ROLL

Instrument for physiotherapy exercises used to improve the muscular flexibility and for relaxing self-massages to relieve muscular tension and favour mental and physical wellbeing.



MODELS	DIMENSIONS	CODE
Cylinder	Ø 15 cm - L 90 cm	1418131
Small cylinder	Ø 15 cm - L 45 cm	1418132



ACCESSORIES

NEEDLE PUMP - Code 1411160

With flexible extension.



Code 1411160

SKIPPING ROPE - Code 1417501
Cotton skipping rope for rhythmic gymnastics. The central section is thicker due to a soft plastic covering. With wooden handles.
Length: 280 cm.


Code 1417501

SKIPPING ROPE - Code 1417502
Coloured cotton skipping rope for rhythmic gymnastics with lead core.
Length: 280 cm.


Code 1417502

MINI TRAMPOLINE - Code 1418111
Mini trampoline.

TECHNICAL SPECIFICATIONS	
Max weight	100 kg
Diameter	98 cm
Height	20 cm


Code 1418621

"BAUMANN" SUPPORTS
Pair of "Baumann" supports, made of clear varnished hardwood.

MODELS	DIMENSIONS	CODE
Medium	168 x 92 x 65 mm	1416081
Small	150 x 88 x 65 mm	1416082



"BAUMANN" BLOCKS - Code 1416083
Pair of "Baumann" blocks, made of clear varnished hardwood.
Dimensions: 140 x 65 x 65 mm.


Code 1416083

MAGIC CIRCLE - Code 1417621
A ring that allows performing several trunk, leg and abdomen exercises. Suitable for arm, back, chest, internal and external thigh and calf toning. With ergonomic handles.
Colour: BLACK Dimensions: Ø 38 cm


Code 1417621





Wooden trolley Code 1418611

WOODEN TROLLEY - Code 1418611
Clear vanished solid wood trolley with three shelves, with castors.
Dimensions: 92 x 68 x 95 cm.



The trolley (code 1418611) is not included in the psychomotor kit.

PSYCHOMOTOR KIT - Code 1418621
A combination of instruments such as sticks, cones, hoops, balls, which qualified staff uses to create psychomotor courses or motor programmes in gyms.

N° 01 wooden stick length 70 cm
N° 01 wooden stick length 80 cm
N° 01 wooden stick length 90 cm
N° 01 pack of PVC sticks length 100 cm
N° 04 cones with three holes and four splits, height 30 cm
N° 02 pairs of small wooden clubs
N° 01 pair of medium supports
N° 01 pair of small supports
N° 01 pair of blocks
N° 02 flat-section nylon hoops 60 cm
N° 02 flat-section nylon hoops 70 cm
N° 02 flat-section nylon hoops 80 cm
N° 04 skipping ropes with handles and thicker central section
N° 04 skipping ropes with lead core, assorted colours
N° 04 420 gr balls for rythmic gymnastics
N° 01 2 kg medicine balls
N° 01 pair ankle/wrist weights 1 kg (2x 0,5 kg)
N° 02 psychomotor training sandbags
N° 02 psychomotor training PVC bags
N° 02 green elastic bands (medium) 120x14,5 cm
N° 01 roll cushion dimensions 35x15 cm
N° 01 roll cushion dimensions 50x25 cm

DIAGNOSTICA
DIAGNOSTICS
DIAGNOSTIC



FISIOTECH

FISIOTECH



PODOSCOPE - Code 1511010

Device for the plantar analysis of feet. Varnished beech structure, crystal footboard and double LED lighting on two sides.

TECHNICAL SPECIFICATIONS	
Foot base dimensions	44 x 39 cm
Max working load	135 kg
Dimensions	50 x 48 x 35 h cm



Code 1511010

STEPS FOR STRESS TEST - Code 1312080

Powder-coated steel structure with wooden adjustable step covered in nonslip material.
It is a useful tool for cardiology stress tests.

TECHNICAL SPECIFICATIONS	
Step heights	30-40-50 cm
Max working load	120 kg
Weight	25 kg
Dimensions	70 x 67 x 129 h cm



Code 1312080

PLUMB BOB - Code 1514010

The plumb bob is useful in measuring deviations in posture and precisely and accurately evaluating whether the patient has a posture of scoliotic nature.



Code 1514010

X-RAY FILM VIEWER

Wall mounted white varnished steel sheet instrument with stainless steel film holder. It has cold fluorescent lamps (6500° K) with unbreakable opaline diffuser.

MODELS	CODE
90 x 60 cm	1512010
70 x 43 cm	1512020



EYE CHART - Code 1512030

Wall mounted white varnished steel sheet. It has cold fluorescent lamps (6500° K) with an unbreakable opaline, screen-printed, acrilic diffuser.

TECHNICAL SPECIFICATIONS	
Distance	3 mt
Dimensions	68 x 31 x thickness 9 cm



Code 1512030

SCALE

Professional instrument with clock-like dial. Aluminium and painted metal structure. Non-toxic, nonslip rubber platform. With nonslip feet.

TECHNICAL SPECIFICATIONS	
Platform dimensions	27 x 40 x 8,5 cm
Dial dimensions	Ø 18 cm
Max working load	150 kg
Division	1000 gr

MODELS	CODE
With column h 90 cm	1513010
Without column	1513020



ACCESSORIES

HEIGHT ROD - Code 1513200

Mechanical telescopic height rod with a strong anodised aluminium rod.
It can be mounted on scales or on the wall.



Code 1513200


ARREDAMENTO
FURNITURE
MOBILIER
FISIOTECH


FISIOTECH

FURNITURE

MEDICAL OFFICE

GLASS DOOR CABINETS

Tempered glass door cabinets with four adjustable shelves and hinged doors, painted steel structure.

TECHNICAL SPECIFICATIONS	
Length	80 / 60 cm
Width	40 cm
Height	180 / 160 cm

MODELS	DIMENSIONS	CODE
Two hinged doors	80 x 40 x 180 cm	1611011
One hinged door	60 x 40 x 160 cm	1611012



DESK - Code 1611101

Desk with three drawers, painted steel structure and plywood board.

TECHNICAL SPECIFICATIONS	
Length	130 cm
Width	80 cm
Height	73 cm



Code 1611101

REVOLVING CHAIR - Code 1611201

Revolving chair with functional aesthetical features, with a high backrest in constant contact with your back.

TECHNICAL SPECIFICATIONS	
Depth	63 cm
Width	54 cm
Height	109 cm
Seat depth	44 cm
Seat width	46 cm
Seat height	55 cm



Code 1611201

UPHOLSTERED CHAIR - Code 1611202

Chair with non-deformable padding, painted steel structure.

TECHNICAL SPECIFICATIONS	
Depth	66 cm
Width	71 cm
Height	75 cm
Seat depth	49 cm
Seat width	50 cm
Seat height	42 cm



Code 1611202

REVOLVING COAT RACK - Code 1611301
Standing, revolving coat rack with umbrella stand and drip tray.



Code 1611301

SCREEN - Code 1611401
Screen made of anodized aluminium tube frame.
Washable, self-extinguishing plastic curtains.
Ø 5 cm castors (central ones with brake).

TECHNICAL SPECIFICATIONS	
Total width	160 cm
Height	170 cm
Section width	55 cm



Code 1611401





REVOLVING STOOL
Stool with five spokes. The height is adjustable with a gas spring. The seat is upholstered in a fire-resistant leatherette.
See colours and materials in the chart on the back cover flap.

TECHNICAL SPECIFICATIONS	
Height adjustment	40-55 cm
Base diameter	60 cm
Seat diameter	35 cm
Max load	120 kg

MODELS	CODE
Black frame with castors	1612011
Black frame with feet	1612012
Chrome-plated frame with castors	1612013
Chrome-plated frame with feet	1612014





UPHOLSTERED BENCH
Powder-coated steel structure, padded seat upholstered in class 1IM fire-resistant, washable leatherette. Maxiumum load 150 Kg.

MODELS	DIMENSIONS	CODE
Small	60 x 30 x 45 cm	1613011
Big	120 x 30 x 45 cm	1613012

FURNITURE

COUCH STAIRS AND SEATING

WOODEN STOOL

Clear varnished hardwood beech stool. Maxiumum load 120 Kg.



MODELS	DIMENSIONS	CODE
Small	18 x 37 x 15 cm	1612201
Big	30 x 50 x 45 cm	1612202



WOODEN CHAIR

Clear-varnished beech wood chairs with or without armrests. Maxiumum load 110 Kg.



MODELS	DIMENSIONS	CODE
Without armrests	51 x 53 x 82 cm	1611203
With armrests	56 x 59 x 82 cm	1611204



MOVING MEDICAL TROLLEY

Multi-purpose medical trolley. It can be fully disassembled. RAL 9001 powder-coated steel tube, RAL 7035 powder-coated steel shelves, complete with Ø 5 cm swivel castors.



TECHNICAL SPECIFICATIONS	
Width	50 / 65 cm
Depth	34 / 45 cm
Height	82 cm

MODELS	DIMENSIONS	CODE
Small, two shelves	50 x 34 x 82 cm	1614010
Small, three shelves	50 x 34 x 82 cm	1614020
Big, two shelves	65 x 45 x 82 cm	1614030
Big, three shelves	65 x 45 x 82 cm	1614040



FURNITURE

COUCH STAIRS AND SEATING

MULTI-PURPOSE TROLLEY

Multi-purpose trolley with welded 18/10 (AISI 304) stainless steel 25-mm tube structure.
With two stainless steel pushing handles.
Complete with four swivel castors Ø125 mm and bumpers.

TECHNICAL SPECIFICATIONS	
Width	82 cm
Depth	52 cm
Height	82 cm



Code 1614101

MODELS	CODE
With two shelves	1614101
With two horizontal drawers	1614111



Code 1614111

COUCH STAIRS

Couch stairs with chrome-plated structure and nonslip steps.

TECHNICAL SPECIFICATIONS	
Step depth	20 cm
Step width	38 cm
First step height	22 cm
Second step height	44 cm
Maximum load	150 kg



Code 1612101

MODELS	DIMENSIONS	CODE
Single step	39 x 28 x 22 cm	1612101
Two steps	39 x 52 x 44 cm	1612102



Code 1612102

CHANGING ROOM BENCHES

Simple changing room bench, varnished steel tube frame and wooden planks.

TECHNICAL SPECIFICATIONS	
Height	45 cm
Depth	35 cm

MODELS	CODE
Width 1 mt	1613021
Width 2 mt	1613022



Code 1613021



Code 1613022

CHANGING ROOM BENCHES

Changing room bench with backrest, varnished steel tube frame and wooden planks.

TECHNICAL SPECIFICATIONS	
Height	85 cm
Depth	40 cm

MODELS	CODE
Width 1 mt	1613023
Width 2 mt	1613024



Code 1613023



Code 1613024

CHANGING ROOM BENCHES

Complete changing room bench, varnished steel tube frame and wooden planks.

TECHNICAL SPECIFICATIONS	
Height	180 cm
Depth	40 cm

MODELS	CODE
Width 1 mt	1613025
Width 2 mt	1613026



Code 1613025



Code 1613026

FURNITURE

CHANGING ROOM FURNITURE

WALL CLOTHES RACK

Wall clothes rack with bag rack, varnished steel structure and wooden planks.




Code 1613032

MODELS	CODE
Width 1 mt with 3 hooks	1613031
Width 2 mt with 6 hooks	1613032

SIMPLE CLOTHES RACK

Wooden wall clothes rack.



Code 1613034

MODELS	CODE
Width 1 mt with 3 hooks	1613033
Width 2 mt with 6 hooks	1613034

FOOTREST

Wooden footrest with rubber feet.




MODELS	CODE
Width 1 mt	1613041
Width 2 mt	1613042

ARTICOLI ORTOPEDICI

ORTHOPAEDIC AIDS

MATÉRIEL ORTHOPÉDIQUE




FISIOTECH

CRUTCHES - Code 1711010

Grey varnished aluminium alloy crutches. Fixed arm support in thermoplastic material. Plastic handle with reflector. Nonslip, stain-resistant rubber stopper.
Package 1 pair.
Height adjustable on 10 levels.

TECHNICAL SPECIFICATIONS	
Min/max height ground to handle	68,5 – 91 cm
Height handle to arm support	24 cm
Weight	0,45 kg
Max load	150 kg



Code 1711010

WALKER WITH FOUR STOPPERS - Code 1712010

Foldable aluminium alloy walker.
With four stoppers.

TECHNICAL SPECIFICATIONS	
Height adjustment	80 – 95 cm
Length	45 cm
Width	55 cm
Weight	2,5 kg
Max load	110 kg



Code 1712010

WALKER WITH TWO STOPPERS AND TWO CASTORS - Code 1712020

Foldable aluminium alloy walker.
With two stoppers and two rigid castors.

TECHNICAL SPECIFICATIONS	
Height adjustment	80 – 95 cm
Length	45 cm
Width	55 cm
Weight	2,5 kg
Max load	110 kg



Code 1712020

WALKER - Code 1712110

Big foldable walker with four swivel castors, of which two with total brake.
Equipped with height-adjustable underarm supports.

TECHNICAL SPECIFICATIONS	
Height adjustment	78 – 95 cm
Underarm support height adjustment	112 – 145 cm
Length	72 cm
Width	55 cm
Castors dimensions	8 cm
Max load	150 kg



Code 1712110

ROLLATOR - Code 1712150

Blue varnished steel foldable Rollator walker. The Rollator walker is equipped with seat, lever brake, basket, tray, and cane holder

TECHNICAL SPECIFICATIONS	
Height adjustment	80 – 98 cm
Length	70 cm
Width	62 cm
Upholstered polyurethane seat	L 43 x h 62 cm
Pivoting front castors dimensions	20 x 5 cm
Rigid rear castors dimensions	20 x 5 cm
Weight	12 kg
Max load	150 kg



Code 1712150


HOME FITNESS
ENTRAÎNEMENT MAISON
FISIOTECH


FISIOTECH

HOME FITNESS

EASY ACCESS EXERCISE BIKE





EASY ACCESS EXERCISE BIKE - Code 1811010
Magnetic home fitness exercise bike.

TECHNICAL SPECIFICATIONS	
Resistance system	Magnetic
Tension levels	8 levels manual tension
Flywheel	10 kg
Display	1 window lcd
Functions	Time / speed / distance / calories / pulse / odometer
Programs	No
Heart Rate Monitoring	Hand pulse
Handlebar	Ergonomic, multiple handles, adjustable
Seat	Graduated compression, horizontally and vertically adjustable
Crank arms	3 pieces
Maximum User Weight	120 kg
Distance seat-pedal	Min. 69 cm – Max 96 cm
Net Weight	28 kg
Gross Weight	31 kg
Set-Up Sizes	105 x 51 x 132 cm (LxWxH)
Package Dimensions	90 x 26 x 64 cm (LxWxH)
Transportation Wheels	Yes
Certifications	CE, ROHS, EN957

HOME FITNESS

ELECTROMAGNETIC EXERCISE BIKE





ELECTROMAGNETIC EXERCISE BIKE - Code 1811020
Electromagnetic home fitness exercise bike.

TECHNICAL SPECIFICATIONS	
Resistance system	Electromagnetic power 8V
Tension levels	16 levels electronically tension
Flywheel	12 kg
Display	1 colour window LCD
Functions	Time / speed / distance / calories / pulse / RPM / program diagram / effort / watt
Programs	Manual – 10 pre set (16 levels)
Heart Rate Monitoring	Hand pulse + Wireless receiver
Heart rate chest belt	Optional
Handlebar	Ergonomic, multiple handles, adjustable
Seat	Graduated compression, horizontally and vertically adjustable
Crank arms	3 pieces
Maximum User Weight	130 kg
Distance seat-pedal	Min. 68 cm – Max 95 cm
Net Weight	29 kg
Gross Weight	33 kg
Set-Up Sizes	103 x 52 x 135 cm (LxWxH)
Package Dimensions	86 x 32 x 64,5 cm (LxWxH)
Transportation Wheels	Yes
Certifications	CE, ROHS, EN957

HOME FITNESS

RECUMBENT EXERCISE BIKE





RECUMBENT EXERCISE BIKE - Code 1811030
Magnetic home fitness exercise bike.

TECHNICAL SPECIFICATIONS	
Resistance system	Magnetic
Tension levels	8 levels manual tension
Flywheel	9 kg
Display	1 window lcd
Functions	Time / speed / distance / calories / pulse
Programs	No
Heart rate monitoring	Hand pulse
Heart rate chest belt	No
Handlebar	Ergonomic, antislip
Seat	High absorption, horizontally adjustable, breathable backrest
Crank arms	3 pcs
Maximum User Weight	135 kg
Net Weight	43 kg
Gross Weight	49 kg
Set-Up Sizes	152 x 60,5 x 125,5 cm (LxWxH)
Package Dimensions	150 x 28 x 69 cm (LxWxH)
Transportation Wheels	Yes
Certifications	CE, ROHS, EN957

HOME FITNESS

TREADMILL



TREADMILL - Code 1812010
Home fitness treadmill.

TECHNICAL SPECIFICATIONS	
Motor	2,5 continuous HP / 3,5 peak HP
Speed	1,0-18 km/h with 0,1 km/h increase - 3 speed quick selection
Incline	Power incline / 18 levels / 3 incline quick selection
Display	1 LCD window 123 x 67 mm, blue backlight
Functions	Time / distance / speed / incline / calories / pulse / steps
Programs	Manual / 25 pre-set with speed and incline variation / 3 customisable / 3 HRC / body fat test / program pause function
Heart rate monitoring	Hand pulse + wireless receiver
Heart rate chest belt	standard
Running deck	Folding, easy unfolding system with hydraulic piston
Running surface/cushioning	46 x 140 cm – 4 air shox + 4 absorbing cushions
Roller size	46 mm front, 42 mm rear
Running belt	Non-slip with double-weave cotton canvas
Handrail buttons	Incline / speed, slip grip
Accessories	Bluetooth speaker + tablet support + 2 bottle holders
Transportation wheels	Yes
All-In-One assembly	No
Set-Up sizes	180 x 75 x 138 cm (LxWxH)
Folded size	113 x 75 x 150 cm (LxWxH)
Package dimensions	76,5 x 33,5 x 174 cm (LxWxH)
Maximum user weight	120 kg
Net weight	63 kg
Gross weight	71 kg
Certifications	CE, ROHS, EN957

STEPPER







STEPPER - Code 1813010
Home fitness stepper.

TECHNICAL SPECIFICATIONS	
Resistance	Indipendent hydraulic pistons - 12 levels
Foldable	Yes
Display	1 LCD window
Functions	Time / step / step/minute / calories
Handlebar	Ergonomic, multiform
Set-Up Sizes	75 x 41.5 x 131 cm (LxWxH)
Folded Size	41.5 x 17 x 98 cm (LxWxH)
Package Dimensions	17 x 41.5 x 98 cm (LxWxH)
Maximum User Weight	100 kg
Net Weight	19 kg
Gross Weight	20 kg

HOME FITNESS

ROWING MACHINE



ROWING MACHINE - Code 1814010
Home fitness rowing machine.

TECHNICAL SPECIFICATIONS	
Resistance system	Magnetic
Resistance	Manual with handle on eight levels
Foldable	Yes
Display	1 LCD window
Functions	Time / strokes / distance / calories / total strokes
Flywheel	6 kg
Set-Up sizes	162 x 55 x 87 cm (LxWxH)
Folded size	87 x 55 x 162 cm (LxWxH)
Package dimensions	142 x 53 x 20,5 cm (LxWxH)
Maximum user weight	100 kg
Net weight	21,5 kg
Gross weight	25 kg

TERMS OF SALE

ORDERS
Orders must be sent by e-mail at contact@fisiotech.com or at info@fisiotech.com. Always specify the article code and the color (if applicable).

OPENING HOURS
From Monday to Friday from 8.30 to 12.30 am and from 1.30 to 5.30 pm.

PRICES
Unless otherwise agreed, the prices here listed are EXW; VAT and other charges are at the buyer's charge. Prices may vary in line with the changes in the cost of raw materials and cost of labor.

MINIMUM OF ORDER
€ 120,00.=. For orders below the minimum of order, you will be charged with € 25,00.=.

SHIPMENT
Shipments are carried out by the shipping company the seller deems the most suitable or by the one chosen by the buyer. The goods are always sent at the buyer's risk, even in case of freight prepaid.

COMPLAINTS
Complaints filed more than eight days since the day of receipt will not be accepted. Returns shall be previously arranged with FISIOTECH and, in any case, shall be carriage free.

RIGHT OF WITHDRAWAL
The right of withdrawal applies to natural persons (consumers) who are acting for purposes unrelated to their commercial activity. Therefore, the right of withdrawal is not applicable for purchases made by retailers, companies or natural persons acting for purposes related to their commercial activity. Only consumers have the right to withdraw within 10 days since the delivery of the goods by returning it to the seller, who will proceed with the refund.

PAYMENTS
In case of delays in the payment, interests will be charged and the agreed terms of payment will be changed.

WARRANTY
Our products' warranty covers only the manufacturing defects for 24 months starting from the date of delivery. The warranty does not cover: consumables and items of normal wear – such as upholstery, power supply, foot pedals/hand controls, power cables and bearings –, breakdowns caused by any impact and/or fall, wrong or improper use, use by non-qualified personnel. The warranty is automatically void if the seal of warranty has been damaged or altered. The warranty is void if the device has been tempered with. Warranty repairs are carried out at our operational headquarters and prior authorization. The products to be repaired under warranty must be shipped to the manufacturer freight prepaid and with the original packaging. A copy of the invoice and of the transport document has to be sent along with a written request in order for the products to be repaired.

LITIGATION APPLICABLE LAW - DISPUTES
The place of jurisdiction for any litigation is the law court of Treviso.

LEAD TIME
The lead time specified in our order confirmation is approximate. Any delay does not entitle the customer to any complaint or refund.

N.B. The seller reserves the right to change features and appearances of the products without any notice. Such modifications will always have as primary purpose to improve the quality/price ratio.

coordinamento
studio grafica design

progetto grafico
fabrizio miniutti

stampa
Combigraf


revisione
01 I22

UPHOLSTERY COLOURS


Beige P1


Light grey G2


Taupe P2


Dark grey G3


Yellow G1


Mustard G4


Orange A1


Red R1


Burgundy B2


Violet V4


Dark blue B3


Blue B1



Dark brown M1


Earth M2


Green V1


Forest green V2


Leaf green V3


Black N1

FISIOTECH

via Toniolo, 44/E-F-G
31028 Vazzola - Treviso - Italy
Tel. +39 0438 47 03 42
Fax +39 0438 47 87 62
info@fisiotech.com
www.fisiotech.com